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PRP-X100

Polymeric anion exchange packings for separation of inorganic and organic anions from 10 to 500 ppm. Easily separate the eight common anions (fluoride through sulfate). Good separation of fluoride from the water dip.

Use with organic solvent from 0 to 100 % for elution of hydrophobic anions or column cleaning.

pH stable from 0 to 14. Use conductivity or UV detection.

Eight column internal diameters from 1.0 to 50.8 mm

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• Bulk material available in 10 g, 100 g and 1000 g quantities.

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
PRP-X110 and PRP-X110S

Polymeric anion exchange packings for separation of inorganic and organic anions from 20 ppb to 20 ppm. Easily separate the anions in EPA method # 300.0 Part A. Good separation of fluoride from the water dip.

Use with organic solvent up to 100% for elution of hydrophobic anions or column cleaning.

pH stable from 0 to 14. PRP-X110 for non-suppressed ion chromatography and PRP-X110S columns for chemically suppressed ion chromatography. Compatible with conductivity or UV detectors.

Four column internal diameters from 1.0 to 4.6 mm

Two column materials: 316 stainless steel and PEEK.

Analytical and semiprep/preparative guard columns.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	Manufacturer
.... E158HM030	PRP-X110	PSDVB*	3 µm	spherical	100 Å	415 m ² /g	Hamilton
.... E158HM070	PRP-X110	PSDVB*	7 µm	spherical	100 Å	415 m ² /g	Hamilton
.... E157HM030	PRP-X110S	PSDVB*	3 µm	spherical	100 Å	415 m ² /g	Hamilton
.... E157HM070	PRP-X110S	PSDVB*	7 µm	spherical	100 Å	415 m ² /g	Hamilton

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· Bulk material available in 10 g, 100 g and 1000 g quantities.

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
Conditions: 1.7 mM Sodium Bicarbonate, 1.8 mM
Sodium Carbonate, 0.1 mM Sodium
Thiocyanate. Isocratic. Ambient.

Flow Rate: 2 ml/min

Injection: 100 µl

Detection: Suppressed Conductivity

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RCX-10

Polymeric anion exchange packing for isocratic or gradient separation of carbohydrates. pH stable from 0 to 14. Isocratic separation of mono and disaccharides. Gradient separation of oligomers up to DP8. Compatible with PAD, conductivity and RI detectors.

Three column internal diameters: 2.1 to 4.6 mm.

Two column materials: 316 stainless steel and PEEK.

Guard columns in stainless steel and PEEK.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	Manufacturer
... E152HM070	RCX-10	PSDVB*	7 µm	spherical	100 Å	-	Hamilton

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

RCX-30

Polymeric anion exchange packing for isocratic or gradient separation of complex carbohydrates. pH stable from 0 to 14. Isocratic or gradient modes. Higher ion exchange capacity support for separation of closely related carbohydrates. Compatible with PAD and RI detectors.

Two column internal diameters: 4.1 to 4.6 mm.

Two column materials: stainless steel and PEEK.

Guard column in PEEK.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	Manufacturer
.... E153HM070	RCX-30	PSDVB*	7 µm	spherical	100 Å	-	Hamilton


* PSDVB = Poly(styrene-divinylbenzene) with Trimethyl-ammonium Exchanger

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
Eurokat

Packing materials in the Eurokat family are sulfonated cross-linked polystyrene copolymers. The higher the cross-linkage, the smaller the pores of the carrier. This support is stable over the entire pH range and is used for special applications. Eurokat is available in 3 different ionic forms. The hydrogen form (Eurokat H) is suitable for the determination of organic acids (Ion Exclusion Chromatography) and mixtures of these with sugars and alcohols (e.g. from fermentation processes) using 0.01 N sulfuric acid as the mobile phase. The calcium and lead forms (Eurokat Ca, Eurokat Pb) are the packing materials of choice for the separation of sugars and sugar alcohols. Sugar alcohols show considerably shorter retention times with Eurokat Ca. Here only purified water as mobile phase is sufficient, making this a particularly environmentally friendly application. All Eurokat columns are characterized by an extended lifetime and robustness.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	Cross Linkage	Manufacturer
.... B340EK100	Eurokat H	Hydrogen	10 µm	spherical	-	-	8 %	KNAUER
.... B340EK200	Eurokat H	Hydrogen	20 µm	spherical	-	-	8 %	KNAUER
.... B360EK100	Eurokat Ca	Calcium	10 µm	spherical	-	-	6 %	KNAUER
.... B360EK200	Eurokat Ca	Calcium	20 µm	spherical	-	-	6 %	KNAUER
.... B350EK100	Eurokat Pb	Lead	10 µm	spherical	-	-	6 %	KNAUER
.... B350EK200	Eurokat Pb	Lead	20 µm	spherical	-	-	6 %	KNAUER

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• Bulk material available in 10 g and 100 g quantities.

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2.1.16 Ion Exchanger polymer based

2.1 Packed HPLC Columns and Bulk Material

CarbEx II

CarbEx is a styrol-divinylbenzol copolymer packing material. This cation exchanger packing material is available with the following counter ions: H and Ca forms.

The packing material exhibits a very good chemical and physical stability and can be used at pressures up to 120 bar and at temperatures up to 90 °C. The eluent may contain up to 70% methanol. The eluent should be well filtered and degassed. The sample should likewise also be filtered before injection. Use an inline filter or a short precolumn (30 mm in length, for example) to protect the column.

Do not stress the column immediately with the full flow rate and corresponding high pressure. Increase the flow rate gradually from 0.1 ml/min up to the working flow rate (for example 0.5-0.8 ml/min for a 250 x 8.0 mm column).


The column should be temperature-controlled when in use. The analysis will run faster and better peak shapes and narrower peaks will result, as well as a lower pressure. The ideal working temperature range is 40-65 °C. Above 50 °C does not produce any significant effect on the separation.

Please note that with the Ca form, very clean double-distilled water must be used.

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
HAMILTON
the measure of excellence

L22

2.1.16 Ion Exchanger polymer based

2.1 Packed HPLC Columns and Bulk Material

PRP-X200

Polymeric cation exchange packing for separation of inorganic cations and organic cations from 20 ppb to 200 ppm. Easily separate mono or divalent cations. pH stable from 0 to 9. Use with organic solvent up to 100% for elution of hydrophobic cations and column cleaning. Compatible with conductivity and UV detectors. Six column internal diameters: 1.0 to 10.0 mm.

Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	Manufacturer
.... E140HM100	PRP-X200	PSDVB*	10 µm	spherical	100 Å	-	Hamilton

* PSDVB = Poly(styrene-divinylbenzene) with Sulfonate Exchanger

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HC-75

Polymeric 7.5% cross-linked soft-gel column for cation, ligand exchange separation of carbohydrates. Water mobile phase. Compatible with RI detectors. Hydrogen form columns for USP L17 applications. Calcium form columns for USP L19 applications. Lead form columns for USP L34 applications.

Three column internal diameters: 4.1 mm, 7.8 mm and 10.0 mm.

One column material: 316 stainless steel.

Three forms: calcium, hydrogen and lead.

Hydrogen & carbonate form guard columns.

Packing Code	Packing	Phase	Particle Size	Pore Size	Surface Area	Manufacturer
.... B148HM090	HC-75	PSDVB*	9 µm	gel-type	-	Hamilton

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

PRP-X300

Polymeric ion exclusion packing for separation of similar alcohols and organic acids. Separate structurally similar alcohols. Separate similar organic acids. Use with organic solvent from 0 to 100% for elution of highly retained compounds.

Five column internal diameters: 1.0 mm to 10.0 mm.

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Packing Code	Packing	Phase	Particle Size	Form	Pore Size	Surface Area	Manufacturer
.... E148HM030	PRP-X300	PSDVB*	3 µm	spherical	100 Å	-	Hamilton
.... E148HM070	PRP-X300	PSDVB*	7 µm	spherical	100 Å	-	Hamilton

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

**Please
ask for additional
Hamilton Phases and Columns**