4. Ventile Inhaltsverzeichnis

4. Ventile

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4. Valves Contents

4. Valves

4.1 Sample Injectors / Switching Valves

- Introduction to Manual Valves / Motor Valves
- Sample Injectors
- "MX" Series Automation Modules
- Tips on Using Injectors
- High Pressure Switching Valves

4.2 Sample Injectors / Switching Valves

- Nano-Scale Sample Injectors / Switching Modules
- "MX" Nano Injektor
- "MX" Nano 6-Port Switching Module
- "MX" Nano 10-Port Switching Module

4.3 Sample Injectors / Switching Valves

- Micro-Scale Sample Injectors / Switching Modules
- Sample Injectors for Microbore Columns
- Model 8125
- Model 7520 Internal Sample Chamber

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
4.1 Sample Injectors / Switching Valves Introduction to Manual Valves

4. Valves

Sample Injectors

How to Choose a Sample Injector

Table I compares the characteristics of manual sample injectors and will help you choose the most suitable model.

Table I. Characteristics of Rheodyne Manual Sample Injectors

Type and Capabilities	Scale	Partial Filling Volumes (Range)	Sample Loop Sizes (Range)	Liquid-Contact Materials	Max MPa ¹	Max T °C	MBB ²	Model ³	
Dual Mode Can load the loop by two methods: 1) partial filling - syringe determines volume without wasting sample 2) complete filling - loop determines volume by overfilling loop	Analytical	1 µl - 2.5 ml	2 µl - 5.0 ml	316 SST, Vespel ⁴ , ceramic, PEEK	48	80	Yes	7725, 7725i	
			5 µl - 5.0 ml		48	80	No	7125	
			1 µl - 5.0 ml	2 µl - 10 ml	PEEK, Tefzel ⁴ , ceramic	34	50	Yes	9725, 9725i
	Micro	0.1 µl - 500 µl	5 µl - 1.0 ml	316 SST, Vespel ceramic, PEEK	48	80	No	8125, 9125	
Preparative	100 µl - 10 ml		2.0 ml - 20 ml	316 SST, PEEK	34	50	Yes	3725(i)-038	
				PEEK	28	50	Yes	3725, 3725i	
Single Mode Can load the loop by one method: complete filling - loop determines volume by overfilling loop	Analytical	Not Applicable	5 µl - 5.0 ml	316 SST, Vespel	48	150	No	7000	
			5 µl - 10 ml	PEEK, Tefzel, ceramic	34	50	No	9010	
	Micro	Not Applicable	0.5 µl - 5.0 µl	316 SST, Vespel	48	150	No	7410	
			0.2 µl - 1.0 µl	316 SST, Vespel	48	80	No	7520	

1. This is the maximum pressure in MPa to which the valve can be adjusted. Some models are shipped from the factory set for lower pressures. 1.0 MPa = 10 bar = 145 psi.

2. MBB™ (Make-Before-Break) is a patented Rheodyne design that provides uninterrupted flow when switching between LOAD and INJECT. MBB™ provides zero sample loss, no cross-contamination and greatly reduces transient pressure shocks.

3. Models with an "i" suffix have a built-in position sensing switch. Models 8125 and 9010 each has a built-in switch.

4. Vespel and Tefzel are trademarks of DuPont.

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4.1 Sample Injectors / Switching Valves Introduction to Manual Valves

4. Valves

Sample Injectors

Types and Capabilities

RHEODYNE models ending in -25 are dual mode injectors. Dual mode injectors can use both the partial-filling and the complete-filling method for loading the sample loop. They are variable volume injectors because they allow the loading of various sample volumes. These dual mode injectors, also called front-loading injectors, have a needle port for loading sample built into the handle. The unique RHEODYNE injection port design allows the tip of the needle to connect directly to the sample loop for no sample loss during loading.

RHEODYNE models ending in -10 are single mode injectors. Single mode injectors use only the complete-filling method to load the sample loop. They are called fixed loop injectors as the sample loop size determines the sample volume. These injectors require a Loop Filler Port accessory, as a needle port is not built into the valve handle. There is not a direct connection between the syringe and the sample loop. Therefore, an excess of sample must be used to overfill the loop filler port and completely fill the sample loop.

RHEODYNE models with an "i" suffix are identical to the models with the same numbers but the "i" designates a built-in position sensing switch. The switch provides the chromatograph with a reproducible start signal to mark the injection time in the data system.

The reproducibility of manual sample injectors depends on operator skill, syringe calibration, and the loading method. Partial-filling method is typically reproducible to 1.0% relative standard deviation (RSD). Complete-filling method is reproducible to 0.1% RSD for loops $\geq 5\mu\text{l}$.

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4.1 Sample Injectors / Switching Valves Introduction to Motor Valves

4. Valves

Sample Injectors

Scale, Sample Volume, and Loop Size

Analytical scale models are for conventional columns with samples from 1.0 µl to 5.0 ml. Micro-scale models are for 1.0 mm and 2.0 mm inner diameter columns. Model 8125 (RH) has a sample range of 0.1 µl to 500 µl, and can be used for both analytical and micro columns. Preparative scale models are for columns with diameters from 1 to 10 cm, and operate at high flow rates with samples from 100 µl to 20 ml.

Liquid Contact Materials

All models have a polymeric rotor seal of Vespel™ (pH 0 to 10 tolerance), Tefzel™, or PEEK (both pH 0 to 14). Stators are 316 stainless steel, titanium, or PEEK. Most models have an inert ceramic stator face assembly.

Make-Before-Break (MBB™)

Models incorporating Rheodyne's patented MBB architecture design provide uninterrupted flow when switching between LOAD and INJECT positions. MBB greatly reduces transient pressure shocks and is beneficial for flow-sensitive detectors, fragile columns, and pumps. Models 7725, 9725, 3725, and "i" versions contain the MBB design.

ChromTRAC™ Mapping

Selected Rheodyne manual valves contain the industry standard ChromTRAC Mapping to color-code your fluid connections. You can identify each port by its colored number, which designate the ChromTRAC color for each system component. Simply coordinate the ChromTRAC colored fittings with the port color.

All valves are ChromTRAC ready. The convenience of ChromTRAC color-coding knobs comes with the RheFlex® Fittings you use for all your connections.

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4. Ventile

Probenaufgabeeventile

4.1 Sample Injectors / Switching Valves Introduction to Motor Valves

4. Valves

Sample Injectors

Table II. Characteristics of Rheodyne "MX" Moduls

Type and Capabilities	Scale		Partial Filling Volumes (Range)	Sample Loop Sizes (Range)	Liquid-Contact Materials	Max MPa ¹	Max T °C	MBB ²	Model
<i>Single Mode</i>	Nano	semi automatic	Not Applicable	10.0 nI			75		MX 7984-000
<i>Dual Mode</i>	Analytical	semi automatic	1 µl - 2.5 ml	2 µl - 5.0 ml	316 SST, PEEK PEEK	35	75	Yes	MX 7925-000
			1 µl - 2.5 ml	2 µl - 5.0 ml		35	75	Yes	MX 9925-000

1. This is the maximum pressure in MPa to which the valve can be adjusted. Some models are shipped from the factory set for lower pressures. 1.0 MPa = 10 bar = 145 psi.

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4.1 Sample Injectors / Switching Valves Introduction to Motor Valves

4. Valves

Sample Injectors

"MX" Series Automation Modules

"MX" Module Specifications

Maximum Pressure	345 bar (5000 psi)
Flow Passages	Nano 0.1 mm (0.004") diameter Analytical 0.3 mm (= .012") diameter
Power Requirements	100-120 VAC, 50-60 Hz
Regulatory Compliance	CE Mark
Remote Control	One line contact closure (open = 1, closed = 2)
Operating Temperature	4° - 75° C
Storage Temperature	-40° - 75° C
Dimensions (H x W x D)	102 mm x 76 mm x 1270 mm (4.5" x 3.0" x 5.0")

A Complete Family of Productivity Enhancing Solutions for Laboratory Automation

- Increase Laboratory Productivity
- Expand Laboratory Capabilities
- Increase Reproducibility
- Improve Analytical Results
- Increase Equipment Reliability
- Save Time and Money
- Make Life Easier

Today's HPLC fluid injection and switching applications require speed, precision, and flexibility. Rheodyne's "MX" Series of productivity enhancing solutions for laboratory automation answer these needs. These completely self-contained, electrically actuated valves are available in a variety of flow paths for nano-, micro-, and analytical-scale applications.

"Manual" or Automatic Actuation

"MX" Modules can be operated by push button, allowing them to function as a "manual" valves, or by contact closure for automated remote control. For multiple automated valve applications, a unique "Snap n' Stack" System allows units to be stacked vertically or connected horizontally to conserve bench space and reduce connection volumes.

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4.1 Sample Injectors / Switching Valves Introduction to Motor Valves

4. Valves

"MX" Series Automation Modules

These small, innovative packages and their simple installation make "MX" Modules the solutions of choice for all high-pressure fluid switching and sample injection applications. Many micro-scale applications may utilize the nano-scale modules.

Nano-Scale Modules

With a mere 10 nl internal sample loop and, more importantly, outstanding dispersion characteristics, Rheodyne's "MX" Nano-Injektor is the ultimate solution for nano- as well as micro-scale applications.

All three Nano-Scale "MX" Modules take advantage of Rheodyne's proprietary biocompatible DuraLife™ III technology for greatly extended duty cycles. "MX" nano-dispersion Six- and Ten-Port Switching Modules are ideal for on-line sample preparation and LC/MS column switching with minimal band broadening.

Analytical-Scale Modules

The "MX" Semi-Automatic Injektor combines the ease of use of Rheodyne's industry standard 7725i with the precision actuation of a electrically actuated valve. With its specially designed needle port, sample is loaded directly into the valve. The push of a button injects your sample. By connecting the injektor to a single contact closure, a sample can be loaded during system equilibration and injected automatically when the system is ready.

Outstanding Performance

"MX" Modules are also available in Analytical-Scale Six-Port and Ten-Port Switching Valves with low dispersion. All "MX" Modules incorporate Rheodyne's time tested miniature "Mighty Valve" - the choice of instrument manufacturers worldwide.

Versatile

"MX" Modules are designed to increase productivity, expand laboratory capabilities, save time, and make life easier.

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Siehe <http://www.rheodyne.com/support/product/technotes/>

4.1 Sample Injectors / Switching Valves Introduction to Motor Valves

4. Valves

"MX" Series Automation Modules


Typical Solutions

- Sample Injection
- High Speed Sample Clean Up
- Two-Column Selection
- Column Backflushing
- Alternating Column Regeneration
- MS Solvent Diversion
- High Speed Sample Enrichment
- Multi-Dimensional Peptide Separation

More Informations about

- Achieving Accuracy and Precision with Rheodyne Manual Sample Injectors
- Pressure Drop of Valves and Tubes from 1 to 1000 ml/min
- Protecting HPLC Columns from Particulate Contamination Using Column Inlet Filters

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4. Valves

Tips on Using Injectors

The Model 7125, 7725 and other front-loading models can use both the partial-filling method and the complete-filling method of loading the loop. Following is a summary of procedures.

1) Steps to inject a sample

- a) With the handle in LOAD, insert the syringe into the needle port until it stops.
- b) Dispense the sample; turn the handle rapidly to INJECT.
- c) Remove the syringe.
- d) Wait until you are ready to run the next sample, then return to LOAD.

2) Flush in the INJECT position, not in LOAD

After turning to INJECT, it is OK to remove the syringe, but keep the handle in this position so the loop is continually flushed with mobile phase.

Manual flushing of the needle port after every injection, to prevent cross contamination, is rarely necessary. A patented direct-connection port design connects the tip of the syringe needle directly to the end of the sample loop; there is no connecting passage that traps sample, which can enter the loop when the next sample is loaded.

It is good practice to flush the needle port after every ten or twenty injections. This keeps it filled with liquid, which baths the syringe needle and dilutes any sample that contaminates this region during movement of the syringe into or out of the port. It also keeps the needle port and vent tube filled with solvent, preventing air from inadvertently entering the loop. To flush, use from 0.1 to 1 ml of mobile phase. Do it while still in the INJECT position, so flow goes directly out vent tube and bypasses the loop which has already been flushed by the pump. Appendix E has more information.

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**Basics
Injection Valves
Switching Valves**

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4. Valves

Tipps on Using Injectors


Figure 1. Schematic of sample flow through mobile phase between tubing walls.


Figure 2. Sample mass (observed peak area) vs. volume of sample dispensed from the syringe, in units of loop volumes, injected onto the column from a Rheodyne dual mode injector such as Model 7725.

3) Flush using the Needle Port Cleaner, not a needle

Use the Needle Port Cleaner (a small Teflon® part without a needle attached), Rheodyne part number 7125-054, attached to a luer tip syringe. This flushes the entire length of the port. A fully inserted needle flushes none of it.

• Fluidic Movement in Tubes

Q: "Why can I load only up to half of the volume of the loop in partial-filling method?"

A: Sample occupies 2 µl of loop for every 1 µl loaded from the syringe. For example, 10 µl of sample spreads out over the entire length of a 20 µl loop. Any more sample loaded will overflow the end of the loop and exit out to waste.

Reproducibility is poor because the volume of sample in the loop is different from the known volume originally loaded by your syringe.

Fluid spreads in a parabolic shape through a tube instead of moving in one plug because the velocity is different at the center of the tube than at the walls. The velocity at the center of the tube is twice the average velocity, and near the wall the velocity is almost zero, creating a parabolic shape. This fluidic movement is called laminar flow (see Figure 1).

In dual mode injectors the sample from the syringe needle loads directly into the sample loop. The sample volume is known since there is no sample waste. The laminar flow phenomenon accounts for the shape of the plot as shown in Figure 2. Note that the plot has three regions:

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Tips on Using Injectors

a) Partial-Filling Region

When the volume dispensed is less than half the loop volume, the curve is linear. Sample has not reached the end of the loop. Within this region, performance depends on the syringe and operator.

b) Nonlinear Region

When the volume dispensed is between half the loop volume and about two loop volumes, the curve is nonlinear. Sample is lost from the loop, so reproducibility is poor. If you dispense a volume equal to the loop size, you are in this region of poor performance.

c) Complete-Filling Region

When the volume of sample dispensed is several loop volumes, the loop contains only pure sample undiluted by residual mobile phase. Within this region, reproducibility is highest.

In the single mode injectors, the sample must pass through a connecting passage before it reaches the sample loop. Since some of the sample dispensed from the syringe remains in the connecting passageway, an unknown amount enters the sample loop. Therefore, single mode injectors achieve high reproducibility only by using the complete-filling method.

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Tips on Using Injectors

• How to Properly Install Fittings

RheFlex® Stainless Steel (and titanium)

A good square-cut, burr-free end on tubing is important for use in your system. The square-cut makes a flat surface for the tube to make contact against the bottom of the valve port. Without this contact, a cavity will form that causes high dispersion and peak distortion such as fronting, tailing, or broadening.

To install the tubing

a) Place the nut (1) and ferrule (2) with threaded portion of the nut and tapered portion of the ferrule onto the tubing (3) (see Figure 1A).

b) Insert the tubing into the valve port (4). Confirm that the tubing is bottomed in the valve port as shown in Figure 1B.

c) While firmly pressing down on the tubing, fingertighten the nut as tight as possible.

d) With the Rheodyne Wrench, designed especially for fittings, tighten 1/4 turn past fingertight. Remove the fitting to confirm that the ferrule is swaged onto the tube.

RheFlex® PEEK Fittings and PEEK Tubing

PEEK tubing is flexible and easy to cut with a sharp blade. Confirm that the ends are clean and straight and that the inner passage is fully open. To install the tubing follow steps a–c in the stainless steel section. Fingertightening of 1/16" PEEK fittings is adequate to make a leak-free connection. The 1/8" PEEK fittings need to be



Figure 1. Cut-away view of stainless steel fittings installation.

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
Figure 2. Cut-away view of PEEK fittings installation.

4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

wrench-tightened. The slotted backside of the ferrule (1) is squeezed down onto the tube (2) by the mating conical surface in the nut (3). (see Figure 2). The nut and ferrule can both be reused many times. Unlike ordinary PEEK fittings, Rheodyne's unique RheFlex PEEK design, specifically the angles and surface contacts between the ferrule and nut, prevents the nut from gripping the ferrule and twisting both the ferrule and the tube during tightening. Otherwise, such twisting stresses the PEEK tubing and lowers the pressure rating of the tubing.

• How to Properly Install Sample Loops

Stainless Steel

Stainless steel sample loops are supplied with fittings that are unswaged onto the tube. It is important that the loop be completely bottomed in the injector port before the ferrule is swaged onto the tube. The depth of the tubing holes may vary slightly from port to port and from valve to valve. A fitting made up in one port may leave a small cavity in another port. The cavity causes high dispersion and peak distortion such as fronting, tailing, or broadening. It is good practice to label loop ends so they will be replaced in the same, respective ports that were used in swaging the ferrules. Hint: swaging ferrules separately on each side, into each respective valve port, makes loop installation easier.

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


Figure 1. Cut-away views (A, B, and C) of stainless steel sample loop installation.

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Basics
Injection Valves
Switching Valves

4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

To install the sample loop

a) Take one end of the loop and place the nut (1) and ferrule (2) onto the tubing (3) with the threaded portion of the nut and tapered portion of the ferrule toward the end (see Figure 1A).

b) Insert the tubing into port (4). Confirm that the tubing is bottomed in the valve port as shown in Figure 1a.

c) While firmly pressing down on the tubing, hand-tighten the nut as tight as possible.

d) With the Rheodyne Wrench, designed especially for fittings, tighten 1/4 turn past fingertight. Remove the loop to confirm that the ferrule is swaged onto the tube.

e) Repeat steps a-d with the other end of the loop while the swaged end remains outside the valve port (see Figure 1B).

f) Reinstall each end of the loop to their respective ports (see Figure 1C).

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RheFlex® PEEK Fittings and PEEK Kapillaren

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
Figure 1. Cut-away view of PEEK sample loop installation.

4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

RheFlex® PEEK Fittings and PEEK Tubing

PEEK loop installation requires steps a–c in the stainless steel section above. Fingertightening of PEEK fittings is adequate to make a leak-free connection. The slotted backside of the ferrule (1) is squeezed down onto the tube (2) by the mating conical surface in the nut (3) (see Figure 2). The nut and ferrule can both be reused many times. Unlike ordinary fittings, the unique RheFlex PEEK design, specifically the angles and surface contacts between the ferrule and nut, prevents the nut from gripping the ferrule and twisting both the ferrule and the tube during tightening. Otherwise, such twisting stresses the PEEK tubing and lowers the pressure rating of the tubing.

The ferrule can slide and reposition itself along the tube when the fitting is reinserted into a port. It is important that the PEEK tubing is completely bottomed in the injector port before the fittings are tightened to avoid leaving an undesired cavity. Both stainless steel and PEEK sample loops are listed in section 4.5.

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4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

• How to Select the Right Rotor Seal

The standard rotor seal in many Rheodyne manual valves is made from a Vespel® blend. This polyimide has low wear and high chemical resistance. Vespel tolerates a pH range of 0 to 10. Solutions more basic than pH 10 dissolve Vespel®, which damages the rotor seal. If you use any solutions above pH 10, Rheodyne recommends a PEEK blend rotor seal. PEEK offers a high chemical resistance and versatility and will tolerate the entire pH range from 0 to 14. Tefzel® blend rotor seals may be appropriate for some applications.


Figure 1. pH range of various rotor seal materials. pH colors are for illustration only.

Genuine Rheodyne rotor seals are matchless in performance and product life. For a quarter of a century, they have exceeded the needs and expectations of chromatographers. The Rheodyne rotor seals are products of rigid manufacturing and quality assurance procedures before they are incorporated into the valves or shipped to the customers. Only genuine Rheodyne parts ensure the continued precision performance of Rheodyne valves. The engineers develop exacting product specifications and designs including the factory-installed rotor seal ring, which optimizes rotor seal efficiency. Rheodyne rotor seals must pass the tougher-than-real-world standards of performance and are made from proprietary-blended polymers, formulated specifically for resistance to repetitive chemical and physical stresses of the entire 0 to 14 pH range. Tested under actual laboratory conditions, Rheodyne rotor seals fully meet the demanding requirements of day-to-day manual instrument use as well as the operating conditions found in today's automated laboratories.

Vespel and Tefzel are registered trademarks of E.I. DuPont.

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
Figure 1A. A square-cut needle stops against the stator face assembly.


Figure 1B. The tip of a pointed needle slips into the stator face, and the tip breaks off as the valve rotates.


Figure 2A. Syringe fitted with needle port cleaner (Part # 7125-054) loading a front-loading valve (Model 7725).


Figure 2B. Loading a loop filler port (Part # 7012).

4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

• How to Use Proper Syringe Needles

With front-loading injectors, it is important to use the correct needle when loading the sample loop. An incorrect needle will damage the valve and can cause poor reproducibility. When the needle is too short, the tip will not reach the needle seal. When the needle is too small in diameter, the seal will not grip tightly enough. Needles with a beveled tip can damage the rotor seal and stator face assembly (see Figure 1B).

The needle should be #22 gauge, and 90° point style (square cut end). Model 3725 requires a #16 gauge needle. Never use a beveled, pointed, or tapered needle. Needle specifications are not critical when using a Loop Filler Port to load the sample loop. However, it is important to tighten the needle port fitting around the needle if using a syringe needle with a slightly smaller diameter than 0.7 mm (0.028"). If the loading method used is complete-filling, a syringe without a needle can be used. A syringe fitted with a Needle Port Cleaner can be used with a front-loading valve (Figure 2A) or with a Loop Filler Port (Figure 2B). Be sure to use proper syringe needle.

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• Sample Loop Loading

Partial-Filling

Use the partial-filling method if you need to conserve sample or if you want to vary sample volume frequently. In partial-filling, the syringe sets the volume injected onto the column. There is no sample waste, and the volume injected onto the column is equal to that dispensed from the syringe. Reproducibility is 1.0% relative standard deviation (RSD). The volume of the sample loaded is limited to half the sample loop volume. For example, the most you can load into a 200 µl sample loop is 100 µl (see Figure 1). This limitation is due to the manner in which fluids move in tubes. Fluidic movement in tubes affects reproducibility.


Figure 1. The sample loop can fill up to half the loaded volume in the partial-filling method.

Complete-Filling

Use the complete-filling method if you have a sufficient amount of sample with which to work, if you do not vary sample volume, or if you need high reproducibility. In complete-filling, the loop sets the volume loaded onto the column. You use excess sample (two to five loop volumes) to replace all the mobile phase in the loop (see Figure 2). Change the loop to vary the sample volume. Reproducibility is typically 0.1% RSD for loop sizes > 5 µl. Accuracy is limited as loop volumes are nominal.


Figure 2. The sample loop is filled in excess in complete-filling method.

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4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

Q: "Which method should I use and which Rheodyne sample injectors use this method?"

A: There are two types of injectors available: dual mode and single mode. Dual mode injectors allow both partial- and complete-filling, whereas single mode injectors allow only complete-filling. See Sample Injectors.

If you are collecting experimental data, sample is scarce, and/or you want to use different sample volumes, a dual mode injector with a large volume sample loop is appropriate. Only dual mode injectors allow the partial-filling method with which you can easily vary your volumes (up to half your sample loop volume) by setting the syringe volume. Once you begin routine analysis and/or you have an abundance of sample, either a dual mode or single mode injector is appropriate. Both types of injectors allow the complete-filling method with which you fill the sample loop in excess. Complete filling maximizes the reproducibility of your results.

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4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

• Avoiding Pressure Transients

Air in the sample loop can cause instantaneous system pressure drop that eventually returns to a normal level. Air causes the pressure to drop when the injector moves from the LOAD to the INJECT position. When large sample loops (>100 µl) are partially loaded, air present in the needle port tube is pushed into the sample loop (see Figure 1). Air can also enter the sample loop from siphoning, which occurs when the vent line is higher than the injection port. In either case, upon injection, the system pressure collapses the air bubble, causing pressure to drop momentarily. A pressure drop in the system caused by air results in changes in retention time, artifact peaks, and affects column performance.


Figure 1. Air (white) present in the needle port tube is pushed by the syringe during loading into the sample loop.


Figure 2. Pathway of the flushing mobile phase using the needle port cleaner when the injector is in INJECT.

Pressure drops can be avoided by removing the air in the needle port tube. Do this by flushing about 1 ml of mobile phase with a luer syringe with needle port cleaner. Keep the needle port tube filled with mobile phase by occasional flushing. Adjust the vent line(s) so that the outlet is at the same horizontal level as the needle port (see Figure 2).

For additional injection troubleshooting, **see:**
<http://www.rheodyne.com/support/product/troubleshooting/index.asp>

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4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

• How to Find and Fix Common Leaks

Leaks cause valuable sample loss. Nobody wants that! The key to the valve holding pressure is the integrity of the sealing surfaces. If there is a scratch on the sealing surface or the needle seal in the rotor seal is damaged, a leak may appear. It is also important to realize what appears to be a leak can instead be a result of siphoning. The following are the three most common situations in which fluid leaks occur.


Figure 1. To reform the needle seal, push the eraser end of a pencil against the needle guide.

1. If fluid leaks out of the needle port only while loading the loop (i.e., while pushing down on the plunger of the syringe), the problem is most likely that the needle seal or the needle port fitting in the loop filler port is not gripping the syringe needle tightly enough to prevent leakage around the syringe needle. Tighten the needle seal grip by pushing down on the needle guide (see Figure 1) or replace the needle port fitting to make a tighter grip on the needle. The tightening reduces the hole diameter of the needle seal and port fitting.

2. If fluid leaks continuously from the needle port or vent lines and/or from the stator-to-stator ring interface, the rotor seal and/or stator face assembly needs to be replaced. Scratches on the rotor seal or cracks in the stator face assembly allow mobile phase to escape and cause cross port leakage.

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<http://www.rheodyne.com/product/troubleshooting/index.asp>

4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

3. If fluid leaks from the needle port and/or vent lines but eventually stops, the cause is most likely siphoning and not a leak. Siphoning occurs if the vent lines are lower or higher than the needle port. Adjust the vent line(s) so that the outlet is at the same horizontal level as the needle port to prevent siphoning (see Figure 2).


Figure 2. Needle port level compared to the level of vent line outlet:

A) Siphoning occurs when the vent line outlet is above the needle port level.

B) Siphoning does not occur if the vent line outlet is the same horizontal level as the needle port.

For other leakage or injection troubleshooting: **see:**

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4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

Tips on Using Injectors

• Burst Pressure of PEEK Tubing

PEEK tubing can be used instead of stainless steel tubing in most applications. PEEK is inert to almost all organic solvents and is biocompatible. Unlike metals, plastics are viscoelastic and, therefore, the yield strengths are not well defined. Many factors affect the burst pressure of PEEK tubing. PEEK tubing will burst at a lower pressure when you:

- Increase the inside diameter.
- Increase the temperature
- Increase the time of exposure.
- Increase the concentration of organic solvents.
- Expose tubing to specific solvents. DMSO, THF, and methylene chloride cause the PEEK tubing to swell. Concentrated nitric acid and sulfuric acid weakens the tubing.

Suggested Maximum Pressure for PEEK Sample Loops.

Tube OD	Tube ID	Sample Loop	Water (1)	1:1 ACN/Water	100% ACN	100% THF	100% IPA
1/16"	0.007"	5 µl	>5500	>5500	4000	2500	>5000
1/16"	0.010"	10, 20 µl	5500	5500	3500	2500	5000
1/16"	0.020"	50, 100, 200 µl	4500	4500	3000	1500	4000
1/16"	0.030"	0.5, 1.0, 2.0, 5.0, 10 ml	3500	3500	3000	1500	4000
1/8"	0.062"	2.0, 5.0 ml	4500	4500	4000	2500	3000
1/8"	0.080"	10, 20 ml	>5000	>5000	2000	1000	2500

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Basics
Injection Valves
Switching Valves

4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

High Pressure Switching Valves

Rheodyne offers high pressure manual switching valves to simplify procedures and improve the speed, resolution and sensitivity of HPLC analysis. The switching valves are available in 316 stainless steel and PEEK, with a choice of 1.6 mm (1/16") or 3.2 mm (1/8") ports.

Column Selection

The six-position switching valves are used in column selection. These valves substitute one column for another without the need to manually disconnect the plumbing. This makes it easy to designate a separate column to each analysis. Designated columns eliminate equilibration delays, reduce interferences, and prolong column life. Turning the valve handle selects the column desired for a particular analysis. The columns switched off-line are automatically sealed at both ends.

Column Switching

The two-position switching valves are used to re-route mobile phase during the chromatographic run without changing separation techniques or to perform sequential separations with different columns and/or mobile phases.

Although the Model 7000 is the most commonly used and versatile switching valve, other models have specific uses such as for three-way or four-way switching patterns.

Many models have flow passages available in both standard bore and large bore, designated with an "L" suffix. L Models use 1/16" fittings and tubing but have larger flow passage diameters than non-L models. L models can accommodate higher flow rates. Large bore tubing can be used when the pressure drop must be limited. Large bore valves have a lower pressure drop than standard bore valves when both valve sizes accommodate the same flow rate.

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
4. Ventile

4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4.Valves

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High Pressure Switching Valves

Figure 1: Pressure drop vs. flow rate for Model 7000 and Model 7000L (large-bore) valves; water at 20°C. Experimental measurements: The flow channel is one stator inlet port, one rotor seal groove, one stator outlet port, and two connecting tubes. Solid squares = (1.0 mm 7000L valve) + (two 1.0 mm x 5.0 cm tubes). Open squares = (0.6 mm 7000 valve) + (two 1.0 mm x 5.0 cm tubes). Cross mark = (0.6 mm 7000 valve) + (two 0.5 mm x 5.0 cm tubes). Solid lines are theoretical values for 10 cm long tubes of 1.0 mm and 0.5 mm ID. Pressure drop is in units of psi. 1.0 MPa = 10 bar = 145 psi.

Table III. Specifications of Rheodyne Manual Switching Valves

Model	Stator Passage Diameter	Factory Set Pressure	Maximum Field Set Pressure	Maximum Temperature
7000, 7030, 7040 (SST & Titanium)	0.6 mm (0.024")	34 MPa (340 bar, 5000 psi)	48 MPa (483 bar, 7000 psi)	150°C*
7060 (SST)	0.4 mm (0.016")	34 MPa (340 bar, 5000 psi)	48 MPa (483 bar, 7000 psi)	80°C
7000L, 7030L, 7040L, 7060L (SST & Titanium)	1.0 mm (0.040")	21 MPa (207 bar, 3000 psi)	34 MPa (340 bar, 5000 psi)	150°C* (7060L: 80°C)
7610-400 (SST)	0.6 mm (0.024")	34 MPa (340 bar, 5000 psi)	34 MPa (340 bar, 5000 psi)	50°C
7610-600 (PEEK)	0.6 mm (0.024")	34 MPa (340 bar, 5000 psi)	34 MPa (340 bar, 5000 psi)	50°C
9010, 9030, 9060 (PEEK)	0.4 mm (0.016")	34 MPa (340 bar, 5000 psi)	34 MPa (340 bar, 5000 psi)	50°C
3000, 3030, 3060 (PEEK)	1.0 mm (0.040")	21 MPa (207 bar, 3000 psi)	28 MPa (276 bar, 4000 psi)	50°C
3000-038, 3030-038, 3060-038 (SST)	1.0 mm (0.040")	28 MPa (276 bar, 4000 psi)	34 MPa (340 bar, 5000 psi)	50°C

SST = stainless steel

* Titanium valves have a maximum temperature of 50°C, due to the Tefzel rotor seal.

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4.1 Sample Injectors / Switching Valves Introduction to Manual Valves / Motor Valves

4. Valves

High Pressure Switching Valves

• Effects of Valves and Tubing on Resolution

The effect of tubing on analytical and micro-scale analyses can be significant. Since dispersion caused by tubing is proportional to the fourth power of diameter, large bore tubing should be avoided when performing analytical scale or micro-scale analyses. A size $\leq 0.25 \text{ mm}$ ($0.010''$) is recommended.



Figure 2: These chromatograms show the loss of resolution caused by the addition of two Model 7060 column selection valves when using connection tubes of two different inside diameters. Conditions for all cases: 4.6 mm x 12.5 cm column, 5 μm C-18 packing, 50% acetonitrile in water, 2.0 ml/min, 21°C, 5.0 μl sample partial-filled into a Model 7125 injector, 10 cm x 0.18 mm (0.007") bore injector outlet tube (to column or valve), 10 cm x 0.18 mm bore detector inlet tube (from column or valve), low dispersion 1.0 cm path UV detector cell, 0.2 sec detector time constant. See text for details.

Consider a system with a Rheodyne injector and column switching valves, and analytical columns with small-bore connecting tubing. The chromatograms in Figure 2, made by using a typical analytical chromatograph, show these effects. Scheme A is the control (injector → column → detector) with no valves in the system. In Schemes B and C, two Model 7060 Six-Position Switching Valves were placed side by side (injector → valve #1 → column → valve #2 → detector).

The injector and detector were connected to these valves by the same tubing used in the control. The extra tubing pieces required to connect the valves to the column were a 10 cm length for valve #1-to-column, and a 35 cm length for column-to-valve #2. The diameters of these tubes are indicated in the Figure caption.

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Position 1 (LOAD) Position 2 (INJECT)
Flow path of "MX" Nano-Scale Sample Injektor

4.2 Sample Injectors / Switching Valves Nano-Scale Sample Injectors / Switching Modules

4. Valves

"MX" Nano Injektor

Rheodyne has created the world's smallest injection valve in a ready-to-use module. The "MX" Nano Injektor has a 10 nl internal sample chamber and outstanding low dispersion characteristics. With its 0.1 mm (0.004") flow passages, the "MX" Nano Injektor can be used for micro-scale applications. Maximum pressure is 345 bar (5.000 psi).

Materials Technology

The "MX" Nano Injektor takes advantage of Rheodyne's proprietary bio-compatible DuraLife™ III technology for greatly extended duty cycles in a typical laboratory automation application.

Actuation

The "MX" Nano Injektor can be used in "manual" or automatic mode. A front panel "REMOTE" button allows the injektor to be actuated by contact closure from a controlling instrument. A bright LCD light indicates when the injektor is in the automatic mode. Simply pressing the button switches the module to "manual" operation. Simply load your sample (Position 1) and press the module's forward position selector button and the module injects your sample (Position 2).

Nano-Scale Sample Injektor

Part Number	Description
MX7984-000	"MX" Nano Injektor, 10 nl internal Loop, Biocompatible DuraLife™ III

Nano Sample Injectors
Switching Valves

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
Position 1 Position 2

Flow path of "MX" 6-Port Switching Valve

4.2 Sample Injectors / Switching Valves Nano-Scale Sample Injectors / Switching Modules

4. Valves

"MX" Nano 6-Port Switching Module

Rheodyne's "MX" Nano Six-Port Switching Module is designed to automate a variety of applications. Like the Nano Injector, this module can be operated by manual push button or, in remote mode, automatically by contact closure with a controlling instrument.

Passages and Pressure

Flow passages are 0.1 mm (0.004") in diameter and the module is rated to 345 bar (5.000 psi).

Extended Lifetime

Rheodyne's "MX" Nano Six-Port Switching Module contains DuraLife™ III wetted surfaces for high duty cycle lifetime. Following are typical applications:

- Two-Column Selection
- Column Backflushing
- High Speed Sample Clean Up and Enrichment
- MS Solvent Diversion

Nano-Scale Six-Port Switching Valve

Part Number	Description
MX7980-000	"MX" Nano-Scale, Two-Position, Six-Port, Biocompatible DuraLife™ III

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
Position 1 (LOAD) Position 2 (INJECT)

Flow path of "MX" Nano-Scale 10-Port Switching Valve

4.2 Sample Injectors / Switching Valves
Nano-Scale Sample Injectors / Switching Modules

4. Valves

"MX" Nano 10-Port Switching Module

Rheodyne's "MX" Nano-Scale Ten-Port Switching Module provides a greater range of automated fluid switching solutions than the six-port module.

Passages and Pressure

This module contains the same 0.1 mm (0.004") flow passages, is rated at 5.000 psi, and is available in stainless steel or biocompatible wetted surfaces.

Applications

Typical applications are:

- Alternating Column Regeneration
- High Speed Sample Enrichment
- Multi-Dimensional Proteomic Peptide Separation

Nano-Scale Ten-Port Switching Valve

Part Number	Description
MX7986-000	"MX" Nano-Scale Two-Position, Ten-Port, Biocompatible DuraLife™ III

**Nano Sample Injectors
 Switching Valves**

Notizen

Notes

A large, empty rectangular area with a light gray background, intended for handwritten notes.

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4.3 Sample Injectors / Switching Valves Micro-Scale Sample Injectors / Switching Modules

4. Valves

Sample Injectors for Microbore Columns

Rheodyne offers three micro-scale injectors that are designed for use with microbore columns (1 mm to 3 mm diameter). The following table outlines their different characteristics.

The method used to load the sample in an injector affects the reproducibility of injectors. Model 8125 can be loaded using partial-filling and complete-filling. Models 7410 and 7520 can only use complete-filling method. Model 7410 uses an accessory injector port (P/N 7012 and 9013).

The precision of all three injectors is typically 1 to 5% relative standard deviation (RSD) depending on the sample volume.

In partial-filling, the syringe sets the volume injected onto the column. There is no sample waste, and the volume injected onto the column is equal to that dispensed from the syringe. The volume of the sample loaded is limited to half the sample loop volume. In complete-filling, the loop sets the volume loaded onto the column. You use excess sample to replace all the mobile phase in the loop.

Comparison of Models 8125, 7410 and 7520

Characteristic	Model 8125	Model 7410	Model 7520
Injection Port Location and Type	Front / Built-In	Rear / Accessory	Front / Built-In
Tubing Connections	0.5 mm (0.020") OD	1/16" OD	1/16" OD
Flow Passage Diameters	0.25 mm (0.010")	0.4 mm (0.016")	0.13 mm (0.005")
Sample Loop Location and Range of Sizes	External: 5 - 1000 µl	Internal: 0.5 - 5 µl	Internal: 0.2 - 1 µl
Complete-Filling Method:	5 µl / 15 µl	5 µl / 30 µl	---
Sample Loop Size (nominal volume)	---	2 µl / 25 µl	---
Volume Required	---	1 µl / 25 µl	1 µl / 7 µl
---	0.5 µl / 40 µl	0.5 µl / 4 µl	
---	---	0.2 µl / 3 µl	
Partial-Filling Method:	0.1 - 500 µl	N/A	N/A
Sample Volumes (volume in syringe)			
Dispersion (variance, µl) [value at specific flow rate]			
5 µl (200 µl/min.)	5.4	6.3	---
(400 µl/min.)	5.8	7.2	---
1 µl (50 µl/min.)	0.8	0.6	1.0
(200 µl/min.)	1.0	0.6	2.0
0.2 µl (50 µl/min.)	0.8	---	0.3
(200 µl/min.)	0.8	---	0.4
Maximum Operating Temperature	80° C	150° C	80° C

Micro Sample Injectors
Switching Valves

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
4.3 Sample Injectors / Switching Valves Micro-Scale Sample Injectors / Switching Modules

4. Valves

Model 8125

Model 8125 is Rheodyne's solution to your micro-scale analyses. Made of 316 stainless steel and designed for 1.0 mm (0.04") and 2.0 mm (0.08") micro-bore columns, the adaptable injector is also compatible with analytical columns (3.0 - 5.0 mm, 0.12 - 0.20"). Model 8125's built-in position sensing switch provides the chromatograph with a reproducible start signal.

This versatile injector allows both partial-filling method (reproducibility of 1.0% RSD) and complete-filling method (reproducibility of 0.1% RSD). To save loading time into the 8125's small flow passages, the largest loop recommended for the complete-filling method is 200 µl. This dual-mode capability varies sample volume desirable for your micro-scale analysis.

Micro-scale 8125 sample loops use 0.5 mm (0.020") OD tubing instead of the conventional 1.6 mm (1/16") OD tubing to provide low-dispersion performance. The versatile 8125 can also accommodate 1.6 mm (1/16") OD tubing. The presence of a mixing cavity between the loop and injector port is less likely when using the smaller size tubing. Cavities may cause high dispersion and peak distortion.

Table A compares the improved resolution using the 8125 to analytical scale injectors, such as the 7725. The improvement is greatest with relatively unretained (low k') peaks.

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
Sample Injector		
Part Number	Description	Stator Material
8125	Dual Mode Micro-Scale Injector	Stainless Steel

4.3 Sample Injectors / Switching Valves Micro-Scale Sample Injectors / Switching Modules

4. Valves

Model 8125

Flow paths of the LOAD and INJECT positions of Model 8125 sample injector for micro-scale analyses.

Flow switching occurs at a flat interface between a polymeric rotor seal and a ceramic stator face assembly in both the stainless steel and PEEK models. You can have confidence in the long seal life of this genuine Rheodyne part combination.

A simple, three-step sample injecting operation involves inserting the syringe into the needle port while in the LOAD position and turning the handle to INJECT. The sample is on its way through your system and when the handle returns to LOAD, the injector is ready for the next injection. The Figure illustrates the flow paths of the positions.

A 5 µl sample loop is furnished.

Specifications

- Max. Operating Pressure: 48 MPa (483 bar, 7.000 psi)
- Max. Operating Temperature: 80° C
- Flow Passage Diameters: 0.3 mm (0.013") and 0.3 mm (0.012")
- Wetted Surfaces: stainless steel and an inert polymer
- Built-in position sensing switch

Table A: Comparison of Observed Column Plates of Rheodyne Analytical and Micro-Scale Injectors

	7725	8125	^
k' = 0.6	2930	5054	72%
k' = 1.5	4653	6904	48%
k' = 7.9	7875	8305	5.0%

UV detector: 1 µl volume, 4 mm path. Sample volume: 2 µl, partial-filling method.
Column: 2 mm ID x 100 mm long, 4 µm C-18. True plates of column = 11.570.

Micro Sample Injectors
Switching Valves

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
4. Ventile

4.3 Sample Injectors / Switching Valves
Micro-Scale Sample Injectors / Switching Modules

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

4.3 Sample Injectors / Switching Valves Micro-Scale Sample Injectors / Switching Modules

4. Valves

Model 7520 Internal Sample Chamber

The 316 stainless steel Model 7520 can inject the smallest volume of all Rheodyne sample injectors with the highest reproducibility. This low-dispersion injector uses a small hole drilled in a flat rotor as an internal sample chamber for precise sample injection. The rotor is available in 0.2, 0.5, and 1 µl. The 0.5 µl rotor is factory installed.

The 7520 requires the complete-filling method of loading the internal sample chamber. Excess sample is required to completely flush mobile phase from the chamber. Reproducibility is 0.1% RSD.

The Figure shows the sliding rotor. The sample loads into the internal sample chamber by inserting a syringe needle into the built-in needle port. The virtually undetectable space between the needle tip and sample chamber contains only a 0.3 µl cavity, which results in little sample waste.

The outlet stator passage is 0.13 mm (0.005") in diameter. A column connection of matching ID is factory installed. The tubing is 1.6 mm (1/16") OD x 5.0 cm (2.0") long. Flow switching occurs at a flat interface between a polymeric rotor seal and a stainless steel stator. You can have confidence in the long seal life of this genuine Rheodyne part combination.

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
Specifications

- Max. Operating Pressure: 48 MPa (483 bar, 7.000 psi)
- Max. Operating Temperature: 80° C
- Flow Passage Diameters: 0.13 mm (0.005")
- Wetted Surfaces: stainless steel and an inert polymer
- Sample Injection Volumes: 0.2 µl, 0.5 µl (standard), 1.0 µl.
Accepts 10-32 male threaded fittings

Flow path of Model 7520
Internal Sample Chamber Micro-Scale Injector.

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Models 7725 / 7725i and 9725 / 9725i

The 316 stainless steel models 7725 and 7725i, and PEEK models 9725 and 9725i are Rheodyne's most advanced manual sample injectors for analytical HPLC.

Specialized features include:

- The Rheodyne patented Make-Before-Break (MBB®) architecture allows continuous flow between LOAD and INJECT positions which greatly reduces transient pressure shocks that disrupts your system. See Figure 1.
- Wide, 30° port angles offer easier access to fittings using the Rheodyne Wrench (Part # 6810).
- Front-end pressure screw makes it easy to adjust and maintain pressure.
- Capability of a reproducible 2 µl sample injection with a 2 µl internal sample loop.
- A built-in position sensing switch ("i" versions) provides the chromatograph with a reproducible start signal.



Figure 1. Flow path of Model 7725 and 9725 with patented Rheodyne MBB design

Figure 1 illustrates the MBB valve design. In the LOAD position, mobile phase flow from pump port to column port travels through both the rotor seal groove and the MBB passage (Position A). As the rotor seal grooves rotate to change from LOAD to INJECT, there is continuous mobile phase flow through both one rotor seal groove and the MBB passage (Position B) until the rotation stops and both rotor seal grooves are connected by the loop. Sample flow begins through the loop to the column just as all flow stops through the MBB passage (Position C). Sample flow never enters the MBB passage.

**Analytical / Preparative
Sample Injectors,
Switching Valves**

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

4.4 Sample Injectors / Switching Valves

Analytical Sample Injectors / Switching Modules

4. Valves

Models 7725 / 7725i and 9725 / 9725i

The injectors contain a patented Rheodyne needle port design that connects the tip of the syringe needle directly to the sample loop ensuring zero sample loss, no cross-contamination, and syringe accuracy. These versatile front-loading injectors allow both partial-filling method (reproducibility of 1.0% RSD) and complete-filling method (reproducibility of 0.1% RSD). This dual mode capability varies sample volumes desirable for your analytical analysis.

Flow switching occurs at a flat interface between a polymeric rotor seal and a ceramic stator face assembly in both the stainless steel and PEEK models. You can have confidence in the long seal life of this

genuine Rheodyne part combination.

A simple, three-step operation involves inserting the syringe into the needle port while in the LOAD position and turning the handle to INJECT. The sample is on its way through your system and when the handle returns to LOAD, the injector is ready for the next injection.

A 20 µl sample loop is furnished. See Table I in 4.1 Introduction to Rheodyne Manual Valves for detailed specifications.

Sample Injectors

Part Number	Description	Stator Material
7725	Dual Mode, Analytical Injector	Stainless Steel
7725i	Dual Mode, Analytical Injector with Switch	Stainless Steel
9725	Dual Mode, Analytical Injector	PEEK
9725i	Dual Mode, Analytical Injector with Switch	PEEK

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
Modelle 7725 / 7725i

Sample Injectors

Part No.	Description, Stator Material
7725	Dual Mode, Analytical Injector, SST
7725i	Dual Mode, Analytical Injector with Switch, SST

Part Number	Description
7010-006	Bearing Ring
7010-015	Isolation Seal
7010-143 (10/pack)	Spring Washers (4)
7125-008	Needle Guide
7125-042	Thrust Bearing
7125-044	Stator Ring
7125-047 (VS)	Rotor Seal
7125-079 (TZ)	Rotor Seal
7125-129 (10/pack)	Handle Set Screws (2)
7125-130	Stator Screws (10/pack)
7725-009	Handle Assembly
7725-010	Stator
7725-026	Stator Face Assembly
7725-034	Pressure Adjusting Nut
8125-005	Body
8125-035	Shaft Assembly
9010-017	60° Stop Ring
9010-018	Position Sensor(7725i only)
9010-090	Stator Locating Pin
9010-091 (10/pack)	Body Locating Pin
9125-123	Needle Port Tube Assembly


4.4 Sample Injectors / Switching Valves Analytical Sample Injectors / Switching Modules

4. Valves

Models 9725 / 9725i / 9125

Sample Injectors

Part No.	Description, Stator Material
9725	Dual Mode, Analytical Injector, PEEK
9725i	Dual Mode, Analytical Injector with Switch, PEEK
9125	Dual Mode, Analytical and Micro Injector, PEEK

Part Number	Description
7010-006	Bearing Ring
7010-015	Isolation Seal
7010-143 (10/pack)	Spring Washers (4)
7125-008	Needle Guide
7125-042	Thrust Bearing
7125-129 (10/pck)	Handle Set Screws (2)
7725-009	Handle Assembly
7725-026	Stator Face Assembly (9725i)
7725-034	Pressure Adjusting Nut
8125-005	Body
8125-035	Shaft Assembly
8125-094	Stator Face Assembly (9125)
9010-016	Stator Ring
9010-017	60° Stop Ring
9010-090	Stator Locating Pin
9010-091 (10/pack)	Body Locating Pin
9125-043 (PK)	Stator
9125-082 (TZ)	Rotor Seal
9125-094	Stator Support Ring
9125-095 (PK)	Rotor Seal
9125-122 (10 pack)	Stator Screws (3)
9125-123	Needle Port Tube Assembly

Analytical / Preparative
Sample Injectors,
Switching Valves

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
Model 7125

Model 7125 is a six-port sample injection valve in which the sample is loaded through a built-in needle port in the front of the valve.

Figure 1 shows the flow diagram of the valve. The circles represent the ports in the valve stator. The dark and white grooves represent the connecting passages in the rotor seal. The large circle represents the needle port.

The sample loop is loaded through the needle port in the LOAD position. Rotation of the knob 60° switches the valve from LOAD to INJECT. In INJECT the mobile phase flows through the loop.


Figure 1. Model 7125 flow diagram.


Figure 2. Guides on adjusting screw and body.


Figure 3. Correct position of rotor seal (slots face the stator).

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

Model 7010 sample loop loading by **pressure loading**. In LOAD, the syringe and Loop Filler Port load the sample loop at port 5 while excess sample exits at port 6. In INJECT (not shown), the valve grooves connect the ports leading from the pump to the column, sending the sample from the sample loop onto the column.

Models 7010 / 7012 Loop Filler Port

Designed specifically for high performance liquid chromatography, the Model 7010 Sample Injection Valve operates up to 7000 psi and injects sample volumes of 5 µl or larger by means of a 1/16" OD sample loop. It provides superior volumetric reproducibility. The Model 7010 valve has six ports and uses the complete-filling method of loading the sample loop..

Model 7012 Loop Filler Port

Model 7012 is an accessory for sample injection valve Models 7010, 7410 and 7413. The Loop Filler Port is a needle insertion port with very small dead volume surrounding the needle when it is inserted.

An adjustable Teflon® seal accommodates needles ranging from 0.5 mm (0.018") to 0.7 mm (0.028") diameter. Conventional microliter syringes with #22 gauge needle can be used. A short length of 0.3 mm (0.012") OD stainless steel tubing connects the needle port to the valve sample loading port.

Specifications

- Max. Operating Pressure: 48 MPa (483 bar, 7.000 psi)
- Max. Operating Temperature: 150° C
- Flow Passage Diameters: 0.6 mm (0.024") and 0.5 mm (0.018")
- Wetted Surfaces: stainless steel and an inert polymer
- Supplied: RheFlex fitting sets for all ports

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
Models 7010 / 7012 Loop Filler Port


Suggested mounting arrangement of Models 7012 and 7010.

Supplied with Loop Filler Port

- Model 7012
 - RheFlex Stainless Steel Fittings set
 - 0.3 mm (0.012") ID Connecting Tube
 - Needle Port Cleaner

- Model 9012
 - RheFlex PEEK Fittings set
 - 0.3 mm (0.010") ID Connecting Tube
 - 0.5 mm (0.020") ID Connecting Tube to rapidly load large sample volumes without high flow resistance
 - Needle Port Cleaner

**Analytical / Preparative
 Sample Injectors,
 Switching Valves**

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

Correct rotor seal orientation as viewed from the stator.

Two-Position Switching Valves
Models 3000-038, 3000, 7610-400,
7610-600 and 9010

Part Number	Description
3000-003	Stator Face Assembly (3000)
3710-008	Rotor Seal (PK) (3000, 3000-038)
3725-006	Stator (3000)
7000-014	Handle Assembly
7010-006	Bearing Ring
7010-015	Isolation Seal
7010-143	Spring Washers (4) (10/pack)
7125-042	Thrust Bearing
7610-017	Stator Face Assembly (7610-400, 7610-600)
7725-034	Pressure Adjusting Nut
8125-005	Body
8125-094	Stator Face Assembly (9010)
9010-016	Stator Ring
9010-017	Stop Ring
9010-051	Rotor Seal (TZ) (9010)
9010-090	Stator Location pin
9010-091	Body Locating Ring
9125-036	Stator (9010)
9125-094	Stator Support Ring (Typ 90 valves only)
9125-122	Stator Screws (Typ 90 valves only, Models 3000, 7610 only)

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
The versatile two-position, six-port and ten-port valves are available in 1/16" and 1/8" port sizes, and 316 stainless steel and PEEK versions. These valves redirect flow among columns during the chromatographic run. They are also useful for selecting between two columns. Compare Model 7000 to the Three-Way Switching Valve (Model 7030) in which each end of the off-line column is independently sealed instead of connected together head-to-tail. Independent seals produce less shock to the column if the valve switches before all the pressure leaves the column. A ten-port valve can often accomplish the same operation that requires two six-port valves.

Model 9010 PEEK sample injector can convert to a six-port switching valve functionally identical to Model 7000 by removing the loop.

Switching Valves

Part Number	Description	Stator Material
7000	Two-Position, Six-Port Switching Valve (1/16")	Stainless Steel
7000L	Two-Position, Six-Port Large Bore Switching Valve (1/16")	Stainless Steel
3000-038	Two-Position, Six-Port Switching Valve (1/8")	Stainless Steel
3000	Two-Position, Six-Port Switching Valve (1/8")	PEEK
7610-400	Two-Position, Ten-Port Switching Valve (1/16")	Stainless Steel
7610-600	Two-Position, Ten-Port Switching Valve (1/16")	PEEK
9010	Single Mode, Analytical Injector	PEEK

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
Three-Way Switching Valves Models 7030(L), 9030, 3030-038 and 3030

Two-position, six-port valves with a double three-way switching pattern are available in 1/16" and 1/8" port sizes, and 316 stainless steel and PEEK versions.

Specifications

- Max. Operating Pressure: 48 MPa (483 bar, 7.000 psi)
L Models - 34 MPa (340 bar, 5.000 psi)
- Max. Operating Temperature: 150° C, Model 7060L - 80° C
- Flow Passage Diameters: Models 7000, 7030, 7040 - 0.6 mm (0.024") and 0.5 mm (0.018");
Models 7060 - 0.4 mm (0.016") and 0.5 mm (0.018");
L Models - 1.0 mm (0.040")
- Wetted Surfaces: stainless steel and an inert polymer

Flow Path of a Three-Way Switching Valve

Switching Valves

Part Number	Description	Stator Material
7030	Three-Way Switching Valve (1/16")	Stainless Steel
7030L	Three-Way Large Bore Switching Valve (1/16")	Stainless Steel
9030	Three-Way Switching Valve (1/16")	PEEK
3030-038	Three-Way Switching Valve (1/8")	Stainless Steel
3030	Three-Way Switching Valve (1/8")	PEEK

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
"MX" Semi-Automatic Analytical-Scale Sample Injector

Easy to Use

The MX Semi-Automatic Injector combines the ease of use of the industry standard 7725i with the precision actuation of a electrically actuated valve.

With its specially designed needle port, sample is loaded directly into the valve. The push of a button injects your sample without the inconsistencies of a manual actuation.

By connecting the injector to a single contact closure, a sample can be loaded during system equilibration and injected automatically when the system is ready.

Flow Path of MX Semi-Automatic Sample Injector

Analytical-Scale Syringe Loading Sample Injector

Part Number	Description
MX7925-000	"MX" Analytical-Scale Syringe Loading Sample Injector, Stainless Steel
MX9925-000	"MX" Analytical-Scale Syringe Loading Sample Injector, Biocompatible

**Analytical / Preparative
 Sample Injectors,
 Switching Valves**

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

Analytical-Scale MX Six-Port Injection Valve

4.4 Sample Injectors / Switching Valves
Analytical Sample Injectors / Switching Modules

4. Valves

"MX" Semi-Automatic Analytical-Scale Sample Injector

Analytical-Scale MX Syringe Loading
Injection Valve

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
Flow Path of "MX" Analytical-Scale Six-Port Switching Valve

4.4 Sample Injectors / Switching Valves
Analytical Sample Injectors / Switching Modules

4. Valves

Analytical-Scale Two-Position, Six-Port Switching Modules

"MX" Analytical Six-Port Switching Module

Rheodyne's "MX" Analytical Six-Port Switching Module is designed to automate a variety of applications. Like the Nano Injektor, this module can be operated by manual push button or, in remote mode, automatically by contact closure by a controlling instrument.

Passages and Pressure

Flow passages are 0.25 mm (0.010") in diameter and the module is rated to 345 bar (5.000 psi).

Extended Lifetime

Following are typical applications:

- Two-Column Selection
- High Speed Sample Clean Up and Enrichment
- Column Backflushing
- MS Solvent Diversion

Analytical-Scale Six-Port Switching Valve

Part Number	Description
MX7900-000	"MX" Two-Position, Six-Port, Stainless Steel
MX9900-000	"MX" Two-Position, Six-Port, Biocompatible

**Analytical / Preparative
 Sample Injectors,
 Switching Valves**

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
Flow Path of "MX" Analytical-Scale Ten-Port Switching Valve

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4. Valves

Analytical-Scale Two-Position, Ten-Port Switching Modules

Rheodyne's "MX" Analytical Ten-Port Switching Module provides a greater range of automated fluid switching solutions than the six-port module.

Passages and Pressure

It contains the same 0.25 mm (0.010") flow passages, is rated at 5.000 psi and takes advantage of Rheodyne's biocompatible DuraLife™ III wetted surfaces.

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For more information see the "MX" Product Bulletin 231 at <http://www.rheodyne.com/PDFs/2321087A.pdf>

For additional illustrated solutions see Rheodyne's "MX" Solutions Guide at <http://www.rheodyne.com/PDFs/2321168A.pdf>

For "MX" Operating manual go to - <http://www.rheodyne.com/PDFs/2321148A.pdf>

Analytical-Scale Ten-Port Switching Valve

Part Number	Description
MX7960-000	"MX" Two-Position, Ten-Port, Stainless Steel
MX9960-000	"MX" Two-Position, Ten-Port, Biocompatible

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
4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

Rotor Seals

The rotor seal is the polymeric disc that makes a high pressure seal against the stator. The seal wears with use and is one of the only parts that may need routine replacement. Stators need replacement only if the ports or sealing surfaces become damaged.

Vespel blend rotor seals have an operating pH range from 0 to 10, Tefzel blend and PEEK blend rotor seals have a pH range from 0 to 14. Strong oxidizing acids such as concentrated nitric and sulfuric are not compatible with PEEK.

Stators are available in 316 stainless steel and PEEK. Our materials of construction have been researched and selected for their physical and mechanical strength.

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4. Ventile

4.5 Accessories / Spare Parts for
Sample Injectors / Switching Valves

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4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

RheBuild® Kits

RheBuild Kits with genuine Rheodyne parts are available for all Rheodyne products. Included in each individualized RheBuild Kit are all parts, tools, and instructions to maintain precision performance of your particular Rheodyne product. RheBuild Kits eliminate individual part ordering.

RheBuild Kits for LabPRO® and EV Automated Fluidics Instrument	
Part Number	Description
1001-999	Model PR100-101
1002-999	Model PR100-102
1005-999	Model PR/EV100-105
1006-999	Model PR/EV100-106
5001-999	Models PR/EV500-101 and PR/EV550-101
5100-999	Models PR/EV500-100 and PR/EV550-100
5104-999	Models PR/EV500-104 and PR/EV550-104
7112-999	Models PR/EV700-112 and PR/EV750-112
7501-999	Models PR/EV700-100 and PR/EV750-100
7502-999	Models PR/EV700-102 and PR/EV750-102
7507-999	Models PR/EV700-107 and PR/EV750-107
7531-999	Models PR703-100 and PR753-100
7004-999	Models PR/EV700-104 and PR/EV750-104

RheBuild Kits for former Models	
Part Number	Description
5700-999	Models 5701, 5703, 5704
7125Ti-999	Model 7125-081

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Needle Port Accessories

Part Number	Description
7012	Stainless Steel Loop Filler Port
9012	PEEK Loop Filler Port
9013	PEEK Needle Port
7125-054	Needle Port Cleaner
9125-076	Suction Needle Adapter (for Model 9725)

Needle Port Accessories

Part Number	Description
7200	2.5 µl Removable Needle Syringe
7200-003	Replacement Needles for 7200, 4/pkg
7215	#22 gauge Stainless Steel Needle for Luer Tip Syringe
3725-056	#16 gauge PEEK Needle for Luer Tip Syringe
3725-086	#16 gauge Stainless Steel Needle for Luer Tip Syringe

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Syringes see chapter 3 section 3.3.

4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

Rheodyne accessories are specially designed for Rheodyne products. Only Rheodyne accessories provide the perfect fit and ensure the precision performance of Rheodyne products.

Needle Port Accessories

Rheodyne's adaptable Loop Filler Ports (Part # 7012 and 9012) are used to load sample from syringe needles or luer tips. The Needle Port (Part # 9013) conserves sample by minimizing the volume between the needle and the valve.

Syringes and Syringe Needles

A wide variety of Syringes and Syringe Needles offers users many sample volumes and needle sizes to precisely load into Rheodyne sample injectors. Some syringes also fit into the Needle Port Cleaner (Part # 7125-054) for flushing out the needle port. All needles have the proper square cut (90° point style) for LC.

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
4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

Mounting Brackets

Rheodyne mounting brackets and panels of different shapes and sizes organize and provide a sturdy support for Rheodyne valves. The Ring Stand Mounting Bracket now allows the valves to mount onto common laboratory equipment.

Position Sensing Switch


The Position Sensing Switch uses a magnetic reed to make a reproducible start signal to record the injection time in the data system of a chromatograph. This convenient, mountable switch is available for sample injectors without "i" versions.

Column Inlet Filters

Column Inlet Filters are useful accessories to extract particles from the sample between the injector and the column. Rheodyne offers Filters and Replacement Filter Discs of different ID sizes for low or medium dispersion. Column Connecting Tubes are also available.

Rheodyne Wrench

The smartly designed Rheodyne Wrench is a double-ended slotted socket wrench that fits over 1/16" and 1/8" ID tubing and loosens and tightens 1/4" and 5/16" stainless steel or PEEK fittings. The "Z" shape of the Rheodyne Wrench provides ideal leverage for changing sample loops and fittings, and keeps one end from restricting the use of the other.

Mounting Bracket Accessories

Part Number	Description
5060	Teflon Rotary Valve Mounting Bracket
7160	Mounting Panel
7160-010	Valve Angle Bracket
7160-029	Ring Stand Mounting Bracket
5060-007	"MX" Ring Stand Mounting Bracket
A4018	Universal Mounting Bracket for mounting manual KNAUER injection valves
A5431	Universal Mounting Bracket Smartline for mounting manual KNAUER injection valve and HPLC column to Smartline instruments

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
4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

Stainless Steel and Titanium Sample Loops

In addition to PEEK sample loops, Rheodyne also supplies a wide variety of 316 stainless steel or titanium sample loops. All sample loops are available with 10-32 threaded nuts and ferrules in 316 stainless steel or titanium.

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4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

Stainless Steel and Titanium Sample Loops

Stainless Steel Loops for 7125 and 7010 Injectors (Do not use for 7725)

Part Number	Description	Bore
7020	5 µl Sample Loop	0.18 mm (0.007")
7021	10 µl Sample Loop	0.30 mm (0.012")
7022	20 µl Sample Loop	0.51 mm (0.020")
7023	50 µl Sample Loop	0.51 mm (0.020")
7024	100 µl Sample Loop	0.51 mm (0.020")
7025	200 µl Sample Loop	0.76 mm (0.030")
7026	500 µl Sample Loop	0.76 mm (0.030")
7027	1.0 ml Sample Loop	0.76 mm (0.030")
7028	2.0 ml Sample Loop	1.0 mm (0.040")
7029	5.0 ml Sample Loop	1.0 mm (0.040")

Stainless Steel Loops for 3725-038 and 3725i-038 Injectors

Part Number	Description	Bore
3065-018	2 ml Sample Loop	2.0 mm (0.080")
3065-019	5 ml Sample Loop	2.0 mm (0.080")
3065-023	10 ml Sample Loop	2.0 mm (0.080")
3065-025	20 ml Sample Loop	2.0 mm (0.080")

Sample loops in other sizes and from other manufacturers on request.

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Part Number	Description
111	Sample Loop, SST, Rheodyne, 5 µl
112	Sample Loop, SST, Rheodyne, 10 µl
113	Sample Loop, SST, Rheodyne, 20 µl
114	Sample Loop, SST, Rheodyne, 50 µl
115	Sample Loop, SST, Rheodyne, 100 µl
116	Sample Loop, SST, Rheodyne, 200 µl
117	Sample Loop, SST, Rheodyne, 500 µl
118	Sample Loop, SST, Rheodyne, 1 ml
119	Sample Loop, SST, Rheodyne, 2 ml
120	Sample Loop, SST, Rheodyne, 3 ml
140	Sample Loop, SST, Rheodyne, 4 ml
141	Sample Loop, SST, Rheodyne, 5 ml

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
4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

Stainless Steel Sample Loops for Rheodyne Valve 7725 and 9725 1/16" OD

Part Number	Description
111-77	Sample Loop, SST, Rheodyne 7725, 5 µl
112-77	Sample Loop, SST, Rheodyne 7725, 10 µl
113-77	Sample Loop, SST, Rheodyne 7725, 20 µl
114-77	Sample Loop, SST, Rheodyne 7725, 50 µl
115-77	Sample Loop, SST, Rheodyne 7725, 100 µl
116-77	Sample Loop, SST, Rheodyne 7725, 200 µl
117-77	Sample Loop, SST, Rheodyne 7725, 500 µl
118-77	Sample Loop, SST, Rheodyne 7725, 1 ml
119-77	Sample Loop, SST, Rheodyne 7725, 2 ml
120-77	Sample Loop, SST, Rheodyne 7725, 3 ml
140-77	Sample Loop, SST, Rheodyne 7725, 4 ml
141-77	Sample Loop, SST, Rheodyne 7725, 5 ml

Material

PEEK or
SST316, DIN1.4401

Tolerances

+/- 0.05 mm (0.002")

Pressure Rating

< 350 bar (< 5000 psi)

Special Information

SST Sample Loops not
recommended for use with
halides

Sample loops in other sizes and from other manufacturers on request.

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4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

PEEK Sample Loops

PEEK sample loops are alternatives to stainless steel loops. The material is inert to almost all organic solvents and is biocompatible which gives PEEK loops added versatility. Rheodyne uses natural PEEK for these sample loops. This means no dye leaching into your sample. Like metal loops, the size designations of PEEK loops are nominal. The actual volumes can differ from the theoretical designations because of the ± 0.05 mm (0.002") tolerance of the tubing bore. Accuracy of large PEEK loops (0.8 mm, 0.030" bore) is about $\pm 14\%$, intermediate loops (0.5 mm, 0.020") $\pm 21\%$, and small loops (0.2 mm, 0.007") $\pm 65\%$.

PEEK Loops for 3725 and 3725i Injectors

Part Number	Description	Bore
3055-018	2 ml Sample Loop	1.6 mm (0.062")
3055-019	5 ml Sample Loop	1.6 mm (0.062")
3055-023	10 ml Sample Loop	2.0 mm (0.080")
3055-025	20 ml Sample Loop	2.0 mm (0.080")

PEEK Loops for 7725, 7725i and PR/EV700-100 Injectors

7123-227	1 μ l Sample Loop (Model PR/EV700-100 only)	internal groove
8021	2 μ l Sample Loop (Models 7725 and 7725i only)	internal groove

PEEK Loops for 9725, 9010, PR/EV750-100 and PR/EV753-100 Injectors Part

Number	Description	Bore
7123-227	1 μ l Sample Loop (Model PR/EV750-100 only)	internal groove
7755-015	2 μ l Sample Loop (9725 only)	internal groove
9055-020	5 μ l Sample Loop	0.18 mm (0.007")
9055-021	10 μ l Sample Loop	0.25 mm (0.010")
9055-022	20 μ l Sample Loop	0.25 mm (0.010")
9055-023	50 μ l Sample Loop	0.51 mm (0.020")
9055-024	100 μ l Sample Loop	0.51 mm (0.020")
9055-025	200 μ l Sample Loop	0.51 mm (0.020")
9055-026	500 μ l Sample Loop	0.76 mm (0.030")
9055-027	1 ml Sample Loop	0.76 mm (0.030")
9055-028	2 ml Sample Loop	0.76 mm (0.030")
9055-029	5 ml Sample Loop	0.76 mm (0.030")
9055-033	10 ml Sample Loop	0.76 mm (0.030")

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Sample loops in other sizes and from other manufacturers on request.

PEEK Physical Strength Characteristics

Although PEEK material is compatible with virtually all solvents, there are many factors that affect burst pressure of PEEK tubing. Factors such as increases in inner diameter, temperature, exposure time, and concentration of organic solvents affect the degradation of PEEK. Other solvents such as THF, methylene chloride, and DMSO cause PEEK tubing to swell while concentrated nitric acid and sulfuric acid weaken tubing.

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Analytical PEEK Sample Loops 1/16" OD

Part Number	Description
2190	Sample Loop, PEEK, Rheodyne, 1 µl
2191	Sample Loop, PEEK, Rheodyne, 2 µl
2192	Sample Loop, PEEK, Rheodyne, 3 µl
2193	Sample Loop, PEEK, Rheodyne, 4 µl
2194	Sample Loop, PEEK, Rheodyne, 5 µl
2196	Sample Loop, PEEK, Rheodyne, 10 µl
2195	Sample Loop, PEEK, Rheodyne, 20 µl
2200	Sample Loop, PEEK, Rheodyne, 30 µl
2201	Sample Loop, PEEK, Rheodyne, 40 µl
2202	Sample Loop, PEEK, Rheodyne, 50 µl
2204	Sample Loop, PEEK, Rheodyne, 100 µl
2206	Sample Loop, PEEK, Rheodyne, 200 µl
2207	Sample Loop, PEEK, Rheodyne, 250 µl
2210	Sample Loop, PEEK, Rheodyne, 500 µl
2222	Sample Loop, PEEK, Rheodyne, 1 ml

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4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

PEEK Sample Loops for Rheodyne Valves

Preparative PEEK Sample Loops 1/8" OD

Part Number	Description
2201-842	Sample Loop, PEEK, Rheodyne, 1 ml
2202-842	Sample Loop, PEEK, Rheodyne, 2 ml
2205-842	Sample Loop, PEEK, Rheodyne, 5 ml
2210-842	Sample Loop, PEEK, Rheodyne, 10 ml
2220-842	Sample Loop, PEEK, Rheodyne, 20 ml
2230-842	Sample Loop, PEEK, Rheodyne, 30 ml
2240-842	Sample Loop, PEEK, Rheodyne, 40 ml
2250-842	Sample Loop, PEEK, Rheodyne, 50 ml

2201-84	Sample Loop, PEEK, Rheodyne, 1 ml
2202-84	Sample Loop, PEEK, Rheodyne, 2 ml
2205-84	Sample Loop, PEEK, Rheodyne, 5 ml
2210-84	Sample Loop, PEEK, Rheodyne, 10 ml
2220-84	Sample Loop, PEEK, Rheodyne, 20 ml
2230-84	Sample Loop, PEEK, Rheodyne, 30 ml
2240-84	Sample Loop, PEEK, Rheodyne, 40 ml
2250-84	Sample Loop, PEEK, Rheodyne, 50 ml

2201-841	Sample Loop, PEEK, Rheodyne, 1 ml
2202-841	Sample Loop, PEEK, Rheodyne, 2 ml
2205-841	Sample Loop, PEEK, Rheodyne, 5 ml
2210-841	Sample Loop, PEEK, Rheodyne, 10 ml

Preparative PEEK Sample Loops 1/16" OD

Part Number	Description
2223-16	Sample Loop, PEEK, Rheodyne, 2 ml
2225-16	Sample Loop, PEEK, Rheodyne, 5 ml
2227-16	Sample Loop, PEEK, Rheodyne, 10 ml
2232-16	Sample Loop, PEEK, Rheodyne, 20 ml

Material

PEEK or
SST316, DIN1.4401

Tolerances

+/- 0.05 mm (0.002")

Pressure Rating

< 350 bar (< 5000 psi)

Special Information

SST Sample Loops not
recommended for use with
halides

Sample loops in other sizes and from other manufacturers on request.

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PEEK and Stainless Steel Sample Loops for Valco Valves

Material

PEEK or SST316, DIN1.4401

Tolerances

+/- 0.05 mm (0.002")

Pressure Rating

< 350 bar (< 5000 psi)

Special Information

SST Sample Loops not recommended for
use with halides

Analytical PEEK Sample Loops 1/16" OD

Part Number	Description
75-5	Sample Loop, PEEK, Valco, 5 µl
75-10	Sample Loop, PEEK, Valco, 10 µl
75-20	Sample Loop, PEEK, Valco, 20 µl
75-50	Sample Loop, PEEK, Valco, 50 µl
75-100	Sample Loop, PEEK, Valco, 100 µl
75-250	Sample Loop, PEEK, Valco, 250 µl
75-500	Sample Loop, PEEK, Valco, 500 µl
75-1000	Sample Loop, PEEK, Valco, 1 ml

Stainless Steel Loops for Valco 1/16" OD

Part Number	Description
75-5	Sample Loop, SST, Valco, 5 µl
75-10	Sample Loop, SST, Valco, 10 µl
75-20	Sample Loop, SST, Valco, 20 µl
75-50	Sample Loop, SST, Valco, 50 µl
75-100	Sample Loop, SST, Valco, 100 µl
75-250	Sample Loop, SST, Valco, 250 µl
75-500	Sample Loop, SST, Valco, 500 µl
75-1000	Sample Loop, SST, Valco, 1 ml
76-500	Sample Loop, SST, Valco, 500 µl
76-1000	Sample Loop, SST, Valco, 1 ml
76-2000	Sample Loop, SST, Valco, 2 ml
76-3000	Sample Loop, SST, Valco, 3 ml
76-4000	Sample Loop, SST, Valco, 4 ml
76-5000	Sample Loop, SST, Valco, 5 ml

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Other dimensions are available on request.

PEEK and SST Sample Loops for Cheminert Valves on request.

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
Abmessungen

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ValvTool

Part Number	Description
800	Tool, ValvTool 1/4" and 5/16" open wrench

Dimensions

for Nuts with 1/4" and 5/16" head

4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

ValvTool

We have designed this unique slotted wrench to tighten nuts where a loop or a capillary may otherwise make it difficult. Not only for Rheodyne and Valco valves and fittings.

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
4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

RheFlex® Stainless Steel Fittings for Cone-Bottom Ports

RheFlex 316 stainless steel fittings sets can be used in Rheodyne stainless steel valves and sample loops. Each 1/16" fittings set contains a 10-32 threaded nut and a ferrule. Three lengths of the 1/16" nut are available: Standard (9.7 mm [0.38"]), Long (13.5 mm [0.53"]), and Extra Long (17.0 mm [0.67"]).

Consult us for titanium fittings.

5/16-24, 1/8" RheFlex® Fittings for Cone-Bottom Ports

RheFlex stainless steel and PEEK fittings are available in 1/8" with nuts containing 5/16-24 threads. These 1/8" fittings sets can be used in all Rheodyne preparative scale valves. Both the 316 stainless steel fittings and RheFlex PEEK fittings must be wrench-tightened.

Due to the unique RheFlex gripping design, the 1/8" ferrule will hold onto PEEK and stainless steel tubing to hold 34 MPa (345 bar, 5.000 psi) without undesirable twisting of the tubing. Replacement 1/8" Stainless Steel and PEEK ferrules are available. A 5/16-24 to 10-32 Adapter is also available. The Adapter allows a 1/16" RheFlex fitting to fit in an 1/8" valve port without causing dead volume.

View the online product product bulletin - <http://www.rheodyne.com/2320626A.html>

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4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

RheFlex® Stainless Steel Fittings and 5/16-24, 1/8" RheFlex® Fittings for Cone-Bottom Ports

RheFlex® 5/16-24, 1/8" Fittings for Cone-Bottom Ports	
Part Number	Description
6000-082	5/16-24, 1/8" RheFlex SST Fittings Set, 1/pkg
6000-083	1/8" RheFlex SST Ferrules, 5/pkg
6000-078	5/16-24, 1/8" RheFlex PEEK Fittings Set, 1/pkg
6000-079	1/8" RheFlex PEEK Ferrules, 5/pkg

Pressure Ratings of RheFlex® Fittings.

RheFlex Fitting	Wrench-tighten MPa (bar, psi)	Finger tighten MPa (bar, psi)
One-Piece Finger Tight	N/A	34 (345, 5000)
1/16" PEEK for Cone-Bottom Ports	48 (483, 7000)	N/A
1/16" Stainless Steel for Cone-Bottom Ports	68 (690, 10,000)	N/A
1/16" Twist-Free for Flat-Bottom Ports	48 (483, 7000) with PEEK tubing	6.8 (68, 1000) with PTFE tubing
1/8" Twist-Free for Flat-Bottom Ports	28 (276, 4000) with PEEK tubing	5.1 (51, 750) with PTFE tubing
1/16" Flangeless for Flat-Bottom Ports	21 (207, 3000) with PEEK tubing	4.7 (47.6, 700) with PTFE tubing
1/8" Flangeless for Flat-Bottom Ports	5.4 (54.4, 800) with PEEK tubing	3.4 (34, 500) with PTFE tubing
1/8" PEEK for Cone-Bottom Ports	34 (345, 5000)	N/A
1/8" Stainless Steel for Cone-Bottom Ports	34 (345, 5000)	N/A

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4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

RheFlex® M4 Fittings for Cone-Bottomed Ports

RheFlex M4 fittings are available in PEEK. Each fitting has a one piece design, which eliminates the need for a separate nut and ferrule. The Rheodyne M4 fitting design provides dependable zero dead volume connections for micro and nano applications. Due to the unique RheFlex gripping design, the M4 fittings will hold onto PEEK or fused silica tubing to hold 34 MPa (345 bar or 5.000 psi). PEEK M4 plugs are also available.

RheFlex® M4 Fittings for Cone-Bottom Ports

Part Number	Description
6000-360	M4 RheFlex PEEK Fittings, 10/pkg
6000-361	M4 RheFlex PEEK Plugs, 10/pkg

RheFlex® Fittings Accessories

Part Number	Description
6000-090	10-32 RheFlex PEEK Plugs, 5/pkg
6510	1/4-28 Coupling, 5/pkg
6518	1/4-28 Plug, 5/pkg
6019	1/4-28 Male to Female Luer Adapter, 1/pkg
6000-076	5/16-24 Male to 10-32 Female Adapter PEEK, 1/pkg
6810	Rheodyne Wrench

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+/- 0,05 mm (0,002")

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< 500 bar (< 7250 psi)

Rheodyne Stainless Steel Nuts and Ferrules

Part Number	Description
56	Nut, SST, short, Rheodyne, 10-32
57	Nut, SST, long head, Rheodyne, 10-32
58	Nut, SST, x-long, Rheodyne, 10-32
59	Ferrule, SST, 1/16"
60	Ferrule, SST, 1/16", 2-step Rheodyne

4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

Rheodyne Type

Our SST Nuts and Ferrules for 1/16" and 1/8" tubing are made of high quality Stainless Steel. We offer nuts (threads 10-32) for 1/16" tubing in three different lengths in order to accommodate any possible connecting situation. These nuts should be used for example with Rheodyne valves.

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
Part Number ZN.5

Part Number ZF.5S6

Valco Stainless Steel Nuts and Ferrules

Part Number	Description
ZN.5	Nut, SST, 1/32"
ZF.5S6	Ferrule, SST, 1/32"
ZN1	Nut, SST, 1/16"
MZN1	Nut, SST, 1/16" medium
LZN1	Nut, SST, 1/16" long
ZF1S6	Ferrule, SST, 1/16"
ZN2	Nut, SST, 1/8"
ZF2S6	Ferrule, SST, 1/8"
ZN4	Nut, SST, 1/4"
ZF4S6	Ferrule, SST, 1/4"

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Pressure Rating

< 500 bar (< 7250 psi)

4.5 Accessories / Spare Parts for Sample Injectors / Switching Valves

4. Valves

Valco Type

Nuts with product numbers with Z are for use with all standard Valco internal fittings and most valves. They may be used with fittings from other manufacturers as well. The L (long) and XL (extra-long) types are for situations where the fitting head may be otherwise inaccessible or where interference between fittings exists, as on many Valco multiposition valves.

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Teflon is a trademark of E.I. DuPont