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
Name	Normal Phase	Reversed Phase	Ionpairing	Application
Silica	x			Non-polar, and moderately polar non-ionic organic compounds.
C1		x		Least retentive of all alkyl group bonded phases for non-polar solutes. Typically used for moderately polar and multi-functional compounds.
C4		x	x	Separation of peptides and proteins. Shorter retention than C8, C18.
C8, Octyl		x	x	Moderately to highly polar compounds, small peptides and proteins, polar pharmaceuticals, steroids, environmental samples.
C18		x	x	Most retentive of the alkyl-bonded phases for non-polar to moderately polar compounds. Used widely for pharmaceuticals, steroids, fatty acids, phthalates, environmental etc.

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1. HPLC in
Laboratory Practice

Phase selection guide

In modern HPLC there is a wide selection of different HPLC column materials available. Each one has been developed to provide a selectivity for a particular type of compound. It can often appear difficult to choose which material will be best for a particular application. The suitability of an HPLC material for an analysis is dependent upon the physical and chemical properties of the analytes and how they will interact with the particular packing.

The following table contains an explanation of the common terms and descriptions in use and an indication of suitable applications for the various chromatographic phases.

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
1.2.1 HPLC Phases and Application areas

1. HPLC in
Laboratory Practice

Name	Normal Phase	Reversed Phase	Ionpairing	Application
CN, Cyano Propyl nitrile	x	x		Unique selectivity for polar compounds, more suitable than base silica for normal-phase gradient separations. When used in reverse-phase the selectivity is different to that of the C8, C18 phases. Useful for a wide range of pharmaceutical applications. Useful for mixtures of very different solutes.
NH ₂ Amino propyl	x	x		RP, Carbohydrate analysis and other polar compounds. Weak anion exchanger, anions and organic acids using buffers and organic modifiers. NP, Alternative selectivity to silica. Good for aromatics.
Phenyl		x	x	Aromatic compounds.
Diol	x	x		RP, Proteins, peptides. NP, similar selectivity to silica, but less polar.
SCX Strong Cation Exchanger				Organic bases.
SAX Strong Anion Exchanger				Organic acids, nucleotides and nucleosides.

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L1

Octadecyl silane chemically bonded to porous silica or ceramic micro-particles, 3 to 10 µm in diameter, or a monolithic rod.

BRAND NAME	MANUFACTURER
Eurosil-Bioselect 300 C18 *	KNAUER
Eurospher 100 C18 *	KNAUER
Hipak ODS AB *	Bischoff Chromatography
Kromasil C18	EKA Nobel
LiChrosorb RP 18	Merck KGaA
LiChrospher 100 RP 18	Merck KGaA
LiChrospher 100 RP 18e	Merck KGaA
Nucleosil C18	Macherey-Nagel
Nucleosil C18 AB	Macherey-Nagel
Partisil ODS	Whatman Inc.
Partisil ODS2	Whatman Inc.
Partisil ODS3	Whatman Inc.
ProntoSIL C18 H	Bischoff Chromatography
ProntoSIL C18 SH	Bischoff Chromatography
ProntoSIL C18 AQ	Bischoff Chromatography
ProntoSIL C18 AQ Plus	Bischoff Chromatography
ProntoSIL C18 ace-EPS	Bischoff Chromatography
ProntoSIL EuroBond C18 *	Bischoff Chromatography

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Octadecyl silane chemically bonded to silica gel of a controlled surface porosity that has been bonded to a solid spherical core, 30 to 50 µm in diameter.

Porous silica particles, 3 to 10 µm in diameter

BRAND NAME	MANUFACTURER
Eurospher 100 Si *	KNAUER
Kromasil Si	EKA Nobel
LiChrosorb Si 100	Merck KGaA
LiChrosorb Si 60	Merck KGaA
LiChrospher Si-60	Merck KGaA

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L4

Silica gel of controlled surface porosity bonded to a solid spherical core, 30 to 50 µm in diameter.

L5

Alumina of controlled surface porosity bonded to a solid spherical core, 30 to 50 µm in diameter.

L6

Strong cation-exchange packing-sulfonated fluorocarbon polymer coated on a solid spherical core, 30 to 50 µm in diameter.

L7

Octylsilane chemically bonded to totally porous silica particles, 3 to 10 µm in diameter.

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L8

An essentially monomolecular layer of aminopropylsilane chemically bonded to totally porous silica gel support, 3 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Eurospher 100 NH2	KNAUER
Kromasil NH2	EKA Nobel
LiChrosorb NH2	Merck KGaA
LiChrospher NH2	Merck KGaA
Nucleosil NH2	Macherey-Nagel
ProntoSIL AMINO	Bischoff Chromatography
ProntoSIL AMINO H	Bischoff Chromatography
ProntoSIL AMINO E	Bischoff Chromatography
Sphere-Image Amino *	Phenomenex
Techsphere Amino *	Grace Vydac
Zorbax NH2	Agilent Technologies

L9

10 µm irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic cation-exchange coating, 3 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Nucleosil SA	Macherey-Nagel
Partisil SCX	Whatman Inc.
Sphere-Image SCX *	Phenomenex
Zorbax SCX	Agilent Technologies

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1.2.2 Corresponding to USP „L“ Phase Listing

L10

Nitrile groups chemically bonded to porous silica particles, 3 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Eurospher 100 CN	KNAUER
Kromasil 60 CN	EKA Nobel
LiChrosorb CN	Merck KGaA
LiChrospher CN	Merck KGaA
Nucleosil CN	Macherey-Nagel
ProntoSIL CN	Bischoff Chromatography
Sphere-Image CN *	Phenomenex
Zorbax CN	Agilent Technologies

L11

Phenyl groups chemically bonded to porous silica particles, 5 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Kromasil Phenyl *	EKA CHEMICALS
Nucleosil Phenyl	Macherey-Nagel
ProntoSIL Phenyl	Bischoff Chromatography
Zorbax Phenyl	Agilent Technologies

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L12

A strong anion-exchange packing made by chemically bonding a quaternary amine to a solid silica spherical core, 30 to 50 µm in diameter.

L13

Trimethylsilane chemically bonded to porous silica particles, 3 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Zorbax TMS	Agilent Technologies

L14

Silica gel having a chemically bonded, strongly basic quaternary ammonium anion-exchange coating, 5 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Nucleosil SB	Macherey-Nagel
Partisil SAX	Whatman Inc.
Sphere-Image SAX *	Phenomenex
Zorbax SAX	Agilent Technologies

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1. HPLC in Laboratory Practice

1.2.2 Corresponding to USP „L“ Phase Listing

Hexylsilane chemically bonded to totally porous silica particles, 3 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
ProntoSIL Spheribond C6	Bischoff Chromatography

Dimethylsilane chemically bonded to porous silica particles, 5 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Eurospher 100 C1 *	KNAUER
Lichrosorb RP2	Merck KGaA
ProntoSIL C1	Bischoff Chromatography

Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the hydrogen form, 7 to 11 µm in diameter.

BRAND NAME	MANUFACTURER
Carbex II (H ⁺) *	VDS optilab
Hamilton HC-75H	Hamilton Co.
Hamilton PRP-X200	Hamilton Co.
Hamilton PRP-X300	Hamilton Co.
Eurokat H *	KNAUER

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L18

Amino and cyano groups chemically bonded to porous silica particles, 3 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Partisil PAC	Whatman Inc..

L19

Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the calcium form, about 9 µm in diameter.

BRAND NAME	MANUFACTURER
Carbex II Ca *	VDS optilab
Hamilton HC 40 *	Hamilton Co.
Hamilton HC 75 Ca	Hamilton Co.
Eurokat Ca *	KNAUER

L20

Dihydroxypropane groups chemically bonded to porous silica particles, 5 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Eurospher 100 Diol *	KNAUER
Kromasil Diol *	EKA CHEMICALS

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L21

A rigid, spherical styrene-divinylbenzene copolymer, 5 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Hamilton PRP-1	Hamilton Co.

L22

A cation-exchange resin made of porous polystyrene gel with sulfonic acid groups, about 10 µm in size.

BRAND NAME	MANUFACTURER
Hamilton PRP-X200	Hamilton Co.
Hamilton PRP-X300	Hamilton Co.

L23

An anion-exchange resin made of porous polymethacrylate or polyacrylate gel with quaternary ammonium groups, about 10 µm in size.

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L24

A semi-rigid hydrophilic gel consisting of vinyl polymers with numerous hydroxyl groups on the matrix surface, 32 to 63 µm in diameter.

L25

Packing having the capacity to separate compounds with a molecular weight range from 100-5000 (as determined by polyethylene oxide), applied to neutral, anionic, and cationic water-soluble polymers. A poly-methacrylate resin base, cross-linked with polyhydroxylated ether (surface contained some residual carboxyl functional groups) was found suitable.

L26

Butyl silane chemically bonded to totally porous silica particles, 3 to 10 µm in diameter.

BRAND NAME	MANUFACTURER
Eurosil Bioselect 300 C4	KNAUER
Eurospher 100 C4 *	KNAUER
Kromasil C4	EKA Nobel
ProntoSIL C4	Bischoff Chromatography

L27

Porous silica particles, 30 to 50 µm in diameter.

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1.2.2 Corresponding to USP „L“ Phase Listing

L28

A multifunctional support, which consists of a high purity, 100 Å, spherical silica substrate that has been bonded with anionic exchanger, amine functionatily in addition to a conventional reversed phase C8 functionality.

L29

Gamma alumina, reverse-phase, low carbon percentage by weight, alumina-based polybutadiene spherical particles, 5 µm in diameter with a pore volume of 80 Å.

L30

Ethyl silane chemically bonded to totally porous silica particles, 3 to 10 µm in diameter.

L31

A hydroxide-selective strong anion-exchange resin-quaternary amine bonded on latex particles attached to a core of 8.5 µm macroporous particles having a pore size of 2000 Å and consisting of ethylvinylbenzene cross-linked with 55% divinylbenzene.

L32

A chiral ligand-exchange resin packing-L-proline copper complex covalently bonded to irregularly shaped silica particles, 5 to 10 µm in diameter.

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L33

Packing having the capacity to separate dextrans by molecular size over a range of 4.000 to 500.000 daltons. It is spherical, silica-based, and processed to provide pH stability.

L34

Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the lead form, about 9 µm in diameter.

BRAND NAME	MANUFACTURER
Hamilton HC-75 Lead Form	Hamilton Co.
Eurokat Pb *	KNAUER

L35

A zirconium-stabilized spherical silica packing with a hydrophilic (diol-type) molecular monolayer bonded phase having a pore size of 150 Å.

L36

A 3.5-dinitrobenzoyl derivatice of L-phenylglycine covalently bonded to 5 µm aminopropyl silica.

L37

Packing having the capacity to separate proteins by molecluar size over a range of 2.000 to 40.000 daltons. It is a polymethacrylate gel.

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1. HPLC in Laboratory Practice

1.2.2 Corresponding to USP „L“ Phase Listing

A methacrylate-based size-exclusion packing for water soluble samples.

L38

A hydrophilic polyhydroxymethacrylate gel of totally porous spherical resin.

L39

Cellulose tris-3.5-dimethylphenylcarbamate coated porous silica particles, 5 µm to 20 µm in diameter.

L40

BRAND NAME	MANUFACTURER
Eurocel 01 *	KNAUER

Immobilized α₁-acid glycoprotein on spherical silica particles, 5 µm in diameter.

L41

Octylsilane and octadecylsilane groups chemically bonded to porous silica particles, 5 µm in diameter.

L42

Pentafluorophenyl groups chemically bonded to silica particles by a propyl spacer, 5 to 10 µm in diameter.

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L44

A multifunctional support, which consists of a high purity, 60 Å, spherical silica substrate that has been bonded with a cationic exchanger, sulfonic acid functionality in addition to a conventional reversed phase C8 functionality.

L45

Beta cyclodextrin bonded to porous silica particles, 5 to 10 µm in diameter.

L46

Polystyrene/divinylbenzene substrate agglomerated with quaternary amine functionalized latex beads about 10 µm in diameter.

L47

High capacity anion-exchange microporous substrate, fully functionalized with a trimethylamine group, 8 µm in diameter.

L48

Sulfonated, cross-linked polystyrene with an outer layer of submicron, porous, anion-exchange micro beads, 15 µm in diameter.

L49

A reversed-phase packing made by coating a thin layer of polybutadiene on to spherical porous zirconia particles, 3 to 10 µm in diameter.

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1. HPLC in Laboratory Practice

1.2.2 Corresponding to USP „L“ Phase Listing

Multifunction resin with reverse-phase retention and strong anion-exchange functionalities. The resin consists of ethylvinylbenzene, 55% cross-linked with divinylbenzene copolymer, 3 to 15 µm in diameter, and a surface area not less than 350 m² per g. Substrate is coated with quaternary ammonium functionalized latex particles consisting of styrene cross-linked with divinylbenzene.

Amylose tris-3,5-dimethylphenylcarbamate-coated, porous, spherical, silica particles, 5 to 10 µm in diameter.

A strong cation exchange resin made of porous silica with sulfopropyl groups, 5 to 10 µm in diameter.

Weak cation-exchange resin consisting of ethylvinylbenzene, 55% cross-linked with divinylbenzene copolymer, 3 to 15 µm diameter. Substrate is surface grafted with carboxylic acid and/or phosphoric acid functionalized monomers. Capacity not less than 500 µEq/column.

A size exclusion medium made of covalent bonding of dextran to highly cross-linked porous agarose beads, about 13 µm in diameter.

A strong cation exchange resin made of porous silica coated with polybutadiene-maleic acid copolymer, about 5 µm in diameter.

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L56

Propyl silane chemically bonded to totally porous silica particles, 3 to 10 µm in diameter.

L62

•
•
•

C30 silane bonded phase on a fully porous spherical silica, 3 to 15 µm in diameter.

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Introduction

A large selection of packing materials have been developed over the years for use in modern HPLC. Just in the area of reversed-phase chromatography (RP-HPLC), more than 700 stationary phases are available. All of these phases were developed to provide optimal selectivity when dealing with specific separation problems. This however also makes it very difficult for the user to find the best material for his separation.

The characteristics of stationary phases are fundamentally dependent on their physicochemical and chemical properties. Important properties include particle diameter, metal content, pore size, surface, the type of functional group on the surface, and also the density of the surface modification (e.g. carbon load). In addition to these parameters, the various stationary phases available can also be differentiated based on their chromatographic properties.

The most important properties of reversed phases are:

- Hydrophobicity
- Silanophilic Activity
- Shape Selectivity
- Polar Selectivity
- Metal Content

These terms will be described briefly below. Each properties' influence on chromatography will be explained through appropriate examples.

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Hydrophobicity

Hydrophobicity refers to the property of a stationary phase to retain a hydrophobic compound as much as possible. The strength of the retention capacity of an individual phase depends substantially on the surface modification and the specific surface of the chromatographic support. C18 phases have the highest hydrophobicity. If the hydrocarbon chains are increased, the loading capacity becomes lower, solely due to steric hinderance. If the hydrocarbon chains are shorter, the analytes will have fewer positions to adsorb, leading to lower hydrophobicity.

The specific surface of the support determines how many functional groups will be bonded. Therefore it is said that the larger the surface, the higher the carbon load, and therefore, the hydrophobicity of the phase.

The surface of a silica gel is dependent on the pore size and pore volume. As long as the pore volume is kept constant, the larger the pores are, the less surface will be available and the lower the hydrophobicity of the stationary phase will be.

This relation is illustrated in Fig. 1 by means of a test mixture separation. The four phases presented are silica gels of different porosities with an identical surface modification.

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
Fig.1 shows the influence of the porosity with the same bonding (C18H)

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Silanophilic Activity

Silanol groups have a large influence on the selectivity of an RP-HPLC phase. On the one hand, undesired secondary interactions result from the silanol groups whereby basic substances will display pronounced tailing through interactions with the acidic silanol groups. On the other hand, such secondary interactions are actually required in some separations to get the desired selectivity. To minimize the silanophilic activity, many stationary phases are "endcapped". This is carried out as an additional treatment step on the already modified (C18, C8, ...) silica surface by treating the silica gel with a smaller silane (usually trisilane). Since this silane is considerably less bulky than the respective C8 or C18 groups, many of the non-reacted silanol groups will react with the silane and therefore become deactivated. Modern columns based on silica gel have very good endcapping and display an exceedingly low silanophilic activity.

Shape Selectivity

The term "shape selectivity" refers to the ability of the stationary phase to differentiate between planar and non-planar molecules (molecular recognition). Highly loaded C18 and C30 phases, in particular, exhibit this property. The effect is based upon the fact that non-planar molecules have a greater space requirement and therefore fewer adsorption positions on the stationary phase as compared to planar analytes. Non-planar molecules are therefore less strongly adsorbed on stationary

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phases with higher shape selectivity. The influence of shape selectivity on the separation of a test mixture is illustrated for three stationary phases in Fig. 2. The shape selectivity of a stationary phase should definitely be taken into consideration when selecting a column.

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Polar Selectivity

The polar selectivity of a stationary phase refers to the property of a phase to undergo polar interactions with analytes. During development of modern RP phases, the newer phases have lost their polar selectivity. Through the synthesis of highly pure silica gels as support materials the surface became more homogeneous, resulting in better coverage by the surface modifications. Modern RP phases with alkyl groups show therefore hardly any polar selectivity. This can be countered however through the targeted introduction of polar functional groups onto the surface of the material or into the alkyl chains.

Fig. 3 compares a chromatogram of a test mixture on a modern C18 reversed phase column with a C18 reversed phase column where an amide group has been introduced into the alkyl chains.

The polar selectivity is of critical importance when selecting a column for method development. If a separation problem cannot be solved on a modern reversed phase material from one manufacturer, it is very probable that the same separation will also not work on another column of the same type from another manufacturer. Therefore the next step should be the use of a column with pronounced polar selectivity.

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Metal Content

Due to the fact that silica gel was earlier produced from sodium silicate, older generations of silica gel possess a very high metal content (predominantly iron, up to 500 ppm). Metals can lead to problems during chromatography of certain analytes. It is known, for instance, that certain proteins cannot be eluted from older HPLC phases. Moreover, the chromatography of complexing agents (e.g. EDTA) with such stationary phases is quite difficult. Since the analytes can only be eluted from the phases as broad asymmetric peaks.

Modern silica gels are produced of tetraethoxysilane and consequently possess metal contents < 10 ppm, making them more suitable for the chromatography of the two examples named above. Fig. 4 shows the separation of bipyridines on two stationary phases of different generations and illustrates the difficulties during chromatography.

Metals in the silica gel also have a direct effect for the acidity of the silica gel. The higher metal content is, the more acidic the silanol groups on the surface of the silica gel will be. This acidity also substantially influences the reactivity of the silica gel. Since the metal content in older generations of silica gel fluctuates and consequently influences the reactivity, the reproducibility of individual batches of stationary phases is affected. Therefore it is recommended to use as much as possible, stationary phases produced from ultra-pure silica gels.

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
1.2 Selection of a Stationary Phase / Column

1. HPLC in Laboratory Practice

1.2.3 Common Criteria to Select the Right Column

Selection of a column according to the properties of the analytes

The selection of the right column for your separation problem depends of course on the properties of your sample. The following properties of an analyte are important for selecting a column: size, functional groups, polarity, spatial structure and matrix.

No special properties

This means: the analyte is either too large or too small (MW between 100 and 5000 g/mol). It possesses neither extremely acidic or extremely basic functional groups. The analytes do not have any significant steric structure and no isomers (e.g. cis/trans isomers) are present. In such a case, a conventional C18 column can be normally used. If possible, a modern C18 phase manufactured from highly pure silica gel should be used. These have, in comparison to older materials, the advantage of better batch to batch reproducibility.

In addition, newer materials also have a better particle size distribution, allowing for higher plate numbers to be achieved at lower back pressure.

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1.2 Selection of a Stationary Phase / Column

1. HPLC in Laboratory Practice

1.2.3 Common Criteria to Select the Right Column

Polar Compounds

In order to be able to retard very polar analytes in reversed phase chromatography, very polar mobile phases are required. Generally it is necessary to work with pure aqueous eluents or even to add ion pair reagents to the mobile phase. Conventional reversed phases cannot be used with pure aqueous eluents. Since the density of the C18 "brushes" on the surface of the silica gel is too high in these phases, the brushes will collapse as soon as the organic portion of the mobile phase is below 5%. In chromatograms, the retention of these analytes fall apart or very broad and asymmetric peaks will be eluted. For aqueous media special stationary phases have been developed. ProntoSIL C18 AQ phases exhibit the same selectivity as conventional C18 phases. ProntoSIL C18 AQ plus phase possesses a polar (hydrophilic) endcapping. The polar groups on the surface attract the polar water and therefore prevent the collapsing of the alkyl chains. ProntoSIL C18 ace-EPS phases are 100% wettable by aqueous eluents. A polar group is attached to the alkyl chains of this phase and therefore produces an analogous behavior as ProntoSIL AQ plus. In comparison to ProntoSIL AQ, the two other possess a polar selectivity imparted by the introduction of a polar functionality, leading to a completely different separation behavior in comparison to conventional C18 phases.

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1.2 Selection of a Stationary Phase / Column

1. HPLC in Laboratory Practice

1.2.3 Common Criteria to Select the Right Column

Non-polar Compounds

Non-polar compounds display a very high retention on conventional reversed phases and can only be eluted from the column with extremely strong eluents. The solution here is to use short-chain reversed phases (C8, C4, C1) which produce substantially weaker interactions with the analytes. Particularly suited to the separation of these compounds are materials with large pores (200 Å and 300 Å) and reversed phases from non-porous silica gel. Since they have less specific surface available they show weaker interactions with the analytes.

An additional possibility is to switch to normal phase mode. Non-polar compounds are often very easy to separate by normal phase chromatography on silica or diol phases. However, you must always deal with retention time shifts in normal phase chromatography since even the slightest change in the water content of the eluents result in an alteration of the retention times. Water is the strongest eluent in normal phase chromatography. Water is also present in all non-polar solvents in trace amounts, making retention time shifts nearly unavoidable.

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1.2 Selection of a Stationary Phase / Column

1. HPLC in Laboratory Practice

1.2.3 Common Criteria to Select the Right Column

High Molecular Weight

As a rule, analytes with high molecular weights have considerably lower diffusion coefficients compared to smaller molecules. Such an analyte's path in and out of the pores is clearly slower in comparison to smaller molecules. This results in broad and asymmetric peaks if the normal pore diameter of 100 Å is used. The use of conventional materials with very large molecular weights can even lead to size exclusion. In such cases, a streak of the compound over the complete separation column is observed in the chromatogram.

For this reason, such compounds are separated on "wide-pore" materials. The typical pore sizes of these supports are about 200 Å and 300 Å. All bondings which are available in smaller pore sizes, are also available in "wide-pore". As an alternative, non-porous materials can also be used. The interaction between the analytes and the stationary phase occur solely on the outer surface of the particle.

The diffusion-controlled migration of the analytes through the pores becomes thereby completely invalidated. Separation columns which have been filled with non-porous support materials offer the most efficient separations for macromolecules.

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1.2 Selection of a Stationary Phase / Column

1. HPLC in Laboratory Practice

1.2.3 Common Criteria to Select the Right Column

Isomers and analytes which differ from each other primarily through their spatial structure

Isomers and analytes which differ from each other primarily through their spatial structure are very difficult to separate. Their hydrophobicity is very similar. Consequently, such separations are regularly carried out by normal phase chromatography on silica phases. The basic problem in normal phase chromatography however is that retention times are shifting. Highly loaded C18 and C30 phases are an alternative and offer an elegant solution to work in reversed phase mode. The distance between the individual alkyl chains on these phases is so minimal, that the separation takes place through the shape selectivity of the phases.

No separation of the analytes on conventional columns

If separation of the analytes is not possible on traditional phases, which is due to the fact that the analytes do not differ enough in terms of polarity and spatial structure, other interactions must be encouraged in order to achieve the desired selectivity. The use of stationary phases with attached polar groups, such as phenyl, cyanopropyl, or aminopropyl groups stimulate the hydrophobic interaction as well as π - π interaction. These modifications achieve completely different selectivities for many separations. Such kind of stationary phases should only be considered as alternatives than as primary choice.

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1.2 Selection of a Stationary Phase / Column

1. HPLC in Laboratory Practice

1.2.3 Common Criteria to Select the Right Column

Special Applications

Special phases have been developed for particular applications in HPLC which possess tailor-made selectivities (see ENVIRO columns). Columns are available for the determination of PAH's, phenols, pesticides, etc. It is recommended that the corresponding column be used for these applications.

A large variety of stationary phases are also available for chiral separations. Since in chiral separations the desired selectivity depends on several parallel interactions. It is very difficult to predict the selectivity in advance. In these cases one should fall back on the applications from the manufacturer or contact the manufacturer for the appropriate selection.

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

1.2.4.1 Narrow-bore and micro-bore HPLC

"Narrow-bore HPLC" involves separation columns with an inner diameter of 2 or 3 mm. If the inner diameter of the column is less than 2 mm, then the column is considered a "micro-bore" column. Recently extremely low volume columns with inner diameters between 50 and 250 μm have been used for HPLC. These columns are called "nano-LC" columns or "capillary HPLC" columns.

Since the volume of the column is proportional to the square of its inner diameter, such low volume HPLC columns can be run with the same linear velocity (i.e. with the same analysis time) with considerably lower flow rates. This results in less eluent consumption and therefore also lowers the cost per analysis. Table 1 compares the eluent consumption of various columns with different inner diameters.

It is clear that by decreasing the column inner diameter, a substantial reduction in eluent consumption can be achieved. The amount of money which is saved by decreasing the amount of high-quality solvents which are necessary, is further increased if the cost of properly disposed solvent waste is considered.

Due to the low volume flow which is given with such columns the coupling to a mass spectrometer is more easy. Given that the eluent must be removed from the sample before it enters the high vacuum of the mass spectrometer, micro-bore columns are particularly practical for such applications.

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1.2 Selection of a Stationary Phase / Column

1. HPLC in
Laboratory Practice

1.2.4 Column Dimensions

Table 1: Eluent consumption and column volume as a function of inner diameter

Innendurchmesser/ Inner diameter [mm]	Flussrate/ Flow rate [μ l/min]	Eluentenverbrauch/ Eluent consumption [ml/h]	Säulenvolumen/ Column volume [μ l]*
4.6	1500	90.0	2740
4.0	1000	60.0	1880
3.0	562.5	33.8	1060
2.0	250.0	15.0	471
1.0	62.5	3.8	118
0.25	3.9	0.2	7
0.1	0.6	0.04	1

* column length: 250 mm, $\epsilon_T = 0.6$

Narrow-and micro-bore columns also have some disadvantages compared to conventional columns. Generally it is not possible to achieve enough precision at low flow rates when using average HPLC pumps. This is a problem particularly with columns of ≤ 1 mm ID. For this reason, flow rates in the lower microliter and nanoliter range are typically achieved through splitting-that is, the eluent is pumped with a normal HPLC pump and then split before the sample injector. The bulk of the eluent flows through a restriction capillary to waste, while a lower constant flow goes to the micro-bore column. In so doing, the advantage of lower eluent consumption is negated.

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

Narrow-bore HPLC columns have the advantage that you can still use these with conventional HPLC equipment without problems, while still achieving a considerable saving in eluent consumption.

However, it is simply not enough to exchange an analytical column (ID 4.6 mm) against a micro-bore column (ID 2.0 mm). Special optimization procedures must be carried out at the HPLC system.

An additional advantage of low-volume columns is their increased mass sensitivity. This means that less sample is required, making analyses possible of very valuable samples. A high mass sensitivity is however not synonymous with higher detection sensitivity. In fact, the opposite is true: Narrow-bore and micro-bore systems are generally less sensitive than conventional analytical systems. Since the column volume has been reduced, the amount injected and the detector cell volume (path-length) has also to be reduced (see chapter on High-Speed HPLC).

A widely accepted thesis is that micro-bore HPLC also gives faster analysis times. This is not the case since by decreasing the column volume, the flow rate has also to be reduced to get the same linear velocity. Due to this fact, the analysis times are identical.

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Schnelle HPLC

1.2 Selection of a Stationary Phase / Column

1. HPLC in
Laboratory Practice

1.2.4 Column Dimensions

1.2.4.2 High-Speed HPLC

In this catalog you will find an assortment of different adsorbents which we offer in various column dimensions. We offer columns with a length of 250, 150 and 125 mm, as well as column inner diameters of 4.6, 4.0, 3.0 and 2.0 mm for analytical HPLC. These columns are usually packed with packing materials in particle sizes of 5 µm or 10 µm.

Shorter columns packed with smaller particles (3 µm or less) are generally used for High-Speed HPLC.

The most important column parameters due to particle sizes and column lengths are shown in the table below.

Table 1:
Characteristic column parameters

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

The values in line 4 of the table shows that the normalized plate count increases significantly if the particle diameter decreases. The column length is the decisive factor in terms of analysis speed. By reducing the column length to 50% the separation efficiency N of the column is also reduced to 50%. By simultaneous reduction of the column length and particle size, the achieved plate count N remains comparable. As shown in line 3 of the table, a typical HPLC separation requires 5.000 - 10.000 plates/column.

With this in mind, you can be sure that short columns packed with small particles are giving the shortest analysis time and the highest separation efficiency. Unfortunately, other system parameters have also to be taken into consideration if column length and particle size are changed. For example, the pressure drop of the HPLC column cannot be chosen unlimited. The back pressure of a column is inverse proportional to the square of the particle diameter (see chapter preparative HPLC).

If a high plate count is required for a particular HPLC separation, a long column packed with larger particles should be used (see line 7 in the table). On the other hand, if analysis speed is of utmost importance, short columns and small particle sizes should be chosen (see line 6 of the table).

An additional advantage of packing materials with smaller particle sizes is that columns packed with such materials can be used at considerably higher flow rates without losing separation efficiency.

Apart from the separation column, the chromatographic process is also influenced by

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Schnelle HPLC

1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

other factors. Band broadening can be caused by diffusion processes in the tubing to the column, as well as in the injector and the detector flow cell (“extra column effects“). It is important to be sure that the extra-column volume of the system does not exceed 20% of the column volume. Since short separation columns are used for High-Speed HPLC, this effect will be particularly effective. Therefore we recommend the use of short columns with inner diameters of 4.0 and 4.6 mm for High-Speed HPLC.

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

1.2.4.3. Preparative HPLC

Introduction

While in analytical HPLC the objective is to identify substances qualitatively and quantitatively, in preparative HPLC the objective is to purify compounds.

In analytical HPLC you work with low concentrations while in preparative HPLC, the column will be overloaded as much as possible due to higher productivity.

Moreover, the objective in preparative HPLC is to dissolve a large amount of compounds in the mobile phase. This allows to use columns packed with larger particles (10 μm or bigger).

You can speak of preparative HPLC if a minimum of 1 g of pure substance can be isolated per run. This can be carried out using columns with an internal diameter > 30 mm. When using columns with an internal diameter from 8 to 30 mm, it is called semi-preparative HPLC.

In HPLC, you distinguish between mass overloading and volume overloading. The loading capacity of an HPLC column depends on the mass of stationary phase, the type of analyte as well as the analyte's solubility in the mobile phase.

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
In Fig. 1 the elution profiles resulting from mass overloading and volume overloading on a preparative column are compared with an analytical separation.

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Fig.1: Comparison of volume overloading to mass overloading

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

Volume overloading

In volume overloading, substances with poor solubility in the mobile phase are injected and elute as a broad zone from the column. Volume overloading has the effect that components elute later in broader elution bands and at lower concentrations. The volume overloading limit is approximately 50 µl/g stationary phase. In practice however, this value is generally exceeded by a factor of 5 to 10. This is valid in Reversed Phase Chromatography. In Normal Phase Chromatography you only will get 10% of the RP-Chromatography values.

Mass overloading

In mass overloading, as much sample as possible is dissolved in the mobile phase and then loaded on the column as concentrated as possible. Since in this case more sample is present than there are adsorption centers available, a part of the sample elutes earlier. The eluting bands will be triangular. The mass overloading limit is about 100 µg/g stationary phase. In practice, this value is typically exceeded by a factor 5 to 10. This is valid in Reversed Phase Chromatography. In Normal Phase Chromatography you only will get 10% of the RP-Chromatography values.

Since the loading capacity for non-charged substances is a 100 times higher than for charged species, acidic and basic substances should be chromatographed as neutral species if possible. This means that preparative HPLC of basic substances should be performed at high pH values and acidic substances at low pH values.

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

Mass overloading should be always preferred over volume overloading in preparative HPLC because the purified products are eluted in higher concentrations and elute earlier from the column. In this case, columns with larger particle diameters can be used without any loss of productivity (see Fig. 2). If a column is overloaded by volume the particle size is of crucial importance. If you use smaller particles you always will get more narrow elution bands.


Fig. 2: Influence of mass overloading on efficiency

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

The influence of particle size in preparative HPLC

In preparative chromatography you should always try to have a cost-effective chromatographic system and to maximize productivity. This is also valid for the separation column. Stationary phases with larger particle diameters are substantially less expensive than those with smaller particle diameters. Larger particles also have a considerably lower pressure drop.

In this case, Darcy's Law is applicable: $\Delta P = \frac{u \cdot \eta \cdot L \cdot \phi}{d_p^2}$

The back pressure of an HPLC column is directly proportional to the linear flow velocity u , the viscosity of the eluent η , the column length L and the specific column resistance ϕ (phi). Inversely proportional to the square of the particle diameter d_p . This means that a column packed with particles size 5 μm has a four-fold back pressure compared to a column packed with 10 μm particles. In addition, you would generally like to increase throughput in preparative HPLC by working with higher flow rates which results in additional back pressure. Due to technical reasons you can only work with a back pressure of up to 15 MPa in preparative HPLC. In analytical HPLC you can work with a back pressure up to 40 MPa or bigger.

Since the column is generally overloaded in preparative HPLC, the separation efficiency of the column plays a smaller role than in analytical HPLC. The optimal particle sizes therefore is in a range of 10 to 15 μm .

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

Upscaling

In chemical production, you normally start with raw material which has to be purified. There are several possibilities. Beside crystallization and distillation, chromatography is gaining importance as the purification method. In case that the more profitable purification methods cannot be utilized, it is possible to use a chromatographic method. In pharmaceutical industry the separation of enantiomeric mixtures by chromatography is one of the best examples.

In those types of situations, the separation has to be first optimized by analytical HPLC. The goal hereby is to choose a packing material and mobile phase allowing the fastest and most selective isocratic separation. To solve your separation problem go for the most selective column. Increasing in selectivity means increase in overloading the column and therefore increase in purified substance per run. Make sure that the stationary phase used is also available in larger particle diameter. This keeps you away to transfer an optimized analytical method to a different packing material, which mostly is difficult.

The analytical column must then be investigated with respect to its loading capacity. The above-mentioned rules of thumb should be applied.

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

The dimension of the preparative column depends on the amount of compound which should be purified per time unit.

The mass (m) or the volume (V) of the compound to purify, as well as the flow rate (F) of the eluent are depending on the ratio of the square of ID (d) from the preparative column (prep) to the square of ID (d) from the analytical column (analytical).

The following equation applies:
$$\frac{F_{\text{analytical}}}{F_{\text{prep}}} = \frac{m_{\text{analytical}}}{m_{\text{prep}}} = \frac{d_{\text{analytical}}^2}{d_{\text{prep}}^2}$$

The characteristics for the most common preparative column dimensions are given in the table below.

F = Flow rate
m = Mass
d = Column inner diameter

Table 1: Upscaling factors for the most common column dimensions.
Linear velocity u ≈ 2 mm/s

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

The linear velocity is kept constant during scale up. Flow rate, injection volume, and mass are modified according to the upscaling factor (see Table 1).

Of course, the instrumental conditions have also to be considered. If scaling up, the internal diameter of the capillaries have to be adapted to the used flow rate (upscaling factor). In contrary to analytical HPLC, the detector should be as insensitive as possible. This is accomplished either by reducing the pass length of the flow cell or by detecting at a wavelength where the compound has the lowest extinction coefficient.

Summary

Ideal process conditions for preparative HPLC :

- Mass overloading
- Solubility of the substances in the mobile phase >20 g/l
- Isocratic conditions (to recycle solvents)
- Particle size of stationary phase 10 to 15 μm
- Potential for high flow rates

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

1.2.4.4 The Column and its Hardware

Besides the selectivity, the performance of the column is of critical importance for the daily work. The most important performance criteria for the quality of a packed column are separation efficiency and stability of the packing. These parameters will depend on the packing material as well as on the column hardware. The unique design of HYPERCHROME column hardware is superior to other hardware solutions in many ways. It's column nuts are perfectly adapted to the column tube, making it possible to use the hardware as a standard column or as a column cartridge. Both types of columns can be used with or without a directly coupled precolumn/pre-column-cartridge. The unique design of the HYPERCHROME HPLC column hardware takes part in the optimization of the chromatographic separation.

The nuts of HYPERCHROME columns may be opened and closed as often as needed. It will not affect the lifetime of the column. Compared to conventional columns, even the packing is protected by a mesh-sandwich at the column ends. Direct distortion or mechanical flow distribution of the packing is thereby eliminated. The mesh-sandwich which replaces typical frits is kept in the column tube by a PTFE seal ring. It also prevents blockages and additional metal contamination which are common for frits. An additional advantage of this column hardware is the unlimited refillability. Quality and performance of a separation will be not influenced by a column tube which has been under Refill. There is no difference in quality between a new column and a refilled column. Before going into the refill process the column tube will be inspected and cleaned like a brand new column tube.

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1.2 Selection of a Stationary Phase / Column

1.2.4 Column Dimensions

1. HPLC in
Laboratory Practice

Reusing column hardware is cost-effective and economical. There is no problem to dispose the used column and you save money compared to buy a new column. In addition you contribute to protect our environment.

Every single column undergoes a quality control test to check its chromatographic performance. This test report is attached to the column so you can be sure that the column meets highest quality criteria.

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Selection Criteria