

5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Flaschen

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Rollrand ND8	Rollrand ND8	Rollrand ND8	Rollrand ND8	Rollrand ND8	Gewinde ND8	Gewinde ND8	Kurzw Gewinde ND9	Gewinde ND10	Rollrand ND11	Rollrand ND11	Rollrand ND11	Rollrand ND11
		<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Screw Neck ND8</i>	<i>Screw Neck ND8</i>	<i>Short Thread ND9</i>	<i>Screw Neck ND10</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>
		LP08 09 0405 LP08 09 0406	LP08 09 0284	LP08 09 0845 LP08 09 0953	LP08 09 0276 LP08 09 0606	LP08 09 0305 LP08 09 0258	LP11 09 0210 LP11 09 0259 LP11 09 0419 LP11 09 0382	LP11 09 0417	LP11 09 0500 LP11 09 0519 LP11 09 0520 LP11 09 0999 LP11 09 0620 LP11 14 1189 LP11 19 0932 LP11 19 1021	LP10 09 0743 LP10 09 1196 LP10 09 1197	LP11 09 0184 LP11 09 0356 LP11 09 0476 LP11 09 0477 LP11 09 0921 LP11 09 0619	LP11 09 0486	LP11 09 0415	LP11 09 0831
Agilent	1050				X		X		X		X			
Agilent	1090				X		X		X					
Agilent	1100						X		X		X			
Agilent	7673A				X				X		X	X		
Agilent	HS7694													
Agilent														
AI	42 vial tray								X		X			X
AI	60 vial tray		X											
AI	CTC A200S	X				X								
AIM							X	X	X		X			X
Alcott	728													
Alcott	738								X					
Alcott									X					X
Antek	736 Unisampler										X			
Antek									X		X			X
Beckman	501				X		X		X		X			
Beckman	502				X		X		X		X			
Beckman	504	X				X			X					
Beckman	507				X		X	X	X		X			
Beckman	502e						X	X	X		X			
Beckman	507e						X	X	X		X			
Beckman	508 (System Gold)								X					
Beckman	Marathon						X		X		X			
Beckman	Promise						X		X		X			
Beckman	Triathlon						X		X		X			
BIO-TEK	HPLC 560 / 565										X			
Bruker	LC51													
CE	42 vial tray				X									
CE	60 vial tray		X		X									
CE	AS 800			X			X	X	X		X		X	X
CE	AS200	X				X								
CE	HS250													
CE	HS500													
CE	HS800													
CE	HS850													
CTC	HS 500													

5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Flaschen

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Schnapping ND11 <i>Snap Ring ND11</i>	Gewinde ND13 <i>Screw Neck ND13</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	
		LP11 09 0627 LP11 09 0644 LP11 09 0645 LP11 14 1190 LP11 19 0933 LP11 19 1022	LP13 09 0222 LP13 09 1335 LP13 09 1336 LP13 09 0280	LP08 14 0641	LP08 14 0513 LP08 14 0595	LP11 14 0544 LP11 14 0545	LP15 14 0548 LP15 14 0562	LP20 09 0342	LP20 09 1405	LP20 09 0297 LP20 09 1223 LP20 09 0440 LP18 09 0864	LP20 09 0795 LP20 09 0796	LP20 09 0873 LP20 09 1222 LP18 09 1306 LP18 09 1310 LP18 09 1307 LP18 09 1311	LP20 09 0801	LP20 09 0802
Agilent	1050													
Agilent	1090													
Agilent	1100													
Agilent	7673A													
Agilent	HS7694										X			
Agilent		X												
AI	42 vial tray	X												
AI	60 vial tray													
AI	CTC A200S													
AIM		X												
Alcott	728			X										
Alcott	738					X								
Alcott														
Antek	736 Unisampler	X												
Antek		X				X								
Beckman	501	X												
Beckman	502	X												
Beckman	504													
Beckman	507	X												
Beckman	502e	X												
Beckman	507e	X												
Beckman	508 (System Gold)								X					
Beckman	Marathon	X												
Beckman	Promise	X												
Beckman	Triathlon	X												
BIO-TEK	HPLC 560 / 565													
Bruker	LC51		X											
CE	42 vial tray													
CE	60 vial tray													
CE	AS 800	X												
CE	AS200													
CE	HS250								X					
CE	HS500								X					
CE	HS800								X					
CE	HS850										X			
CTC	HS 500								X					



5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Vials

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Rollrand ND8	Rollrand ND8	Rollrand ND8	Rollrand ND8	Rollrand ND8	Gewinde ND8	Gewinde ND8	Kurzgewinde ND9	Gewinde ND10	Rollrand ND11	Rollrand ND11	Rollrand ND11	Rollrand ND11	
		<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Screw Neck ND8</i>	<i>Screw Neck ND8</i>	<i>Short Thread ND9</i>	<i>Srew Neck ND10</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	
		LP08 09 0405 LP08 09 0406	LP08 09 0284	LP08 09 0845 LP08 09 0953	LP08 09 0276 LP08 09 0606	LP08 09 0305 LP08 09 0258	LP11 09 0210 LP11 09 0259 LP11 09 0419 LP11 09 0382	LP11 09 0417	LP11 09 0500 LP11 09 0519 LP11 09 0520 LP11 09 0999 LP11 09 0620 LP11 14 1189 LP11 19 0932 LP11 19 1021	LP10 09 0743 LP10 09 1196 LP10 09 1197	LP11 09 0184 LP11 09 0356 LP11 09 0476 LP11 09 0477 LP11 09 0921 LP11 09 0619	LP11 09 0486	LP11 09 0415	LP11 09 0831	
CTC	PAL	X			X	X									
CTC	A 200 S	X				X					X				
Dani	ALS 39.80								X		X				
Dani	ALS 86.80								X		X				
Dani	HS39.50														
Dani	HS86.50														
Dionex	Autoselect AS50						X		X	X	X				
D-Star	DAS-10						X								
Dynatech	42 vial tray								X		X			X	
Dynatech	60 vial tray		X												
Dynatech	LC2000			X											
Finnigan	A 200 S	X				X									
Fisons	60 vial tray		X												
Fisons	AS 800			X			X	X	X		X		X	X	
Fisons	HS500														
Fisons	HS800														
Fisons	HS850														
GBC									X		X				
Gerstel	MPS2	X			X	X									
Gilson	(not ASPEC)	X				X									
Gilson	201/202														
Gilson	221/222														
Gilson	231/401						X		X						
Gilson	232/402						X		X						
Gilson	ASPEC						X		X						
Gilson	HT300A				X			X	X		X				
HTA	HT300A							X	X		X				
ICI	LC1600	X				X			X						X
Infochroma									X		X				
KNAUER	Smartline 3800						X		X		X				
KNAUER	Smartline 3900								X		X				
Kontron	MSI 660								X		X				
Kontron					X			X	X		X		X		
LDC	713-60		X						X		X				
LDC	Marathon						X		X		X				

5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Flaschen

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Schnapping ND11 <i>Snap Ring ND11</i>	Gewinde ND13 <i>Screw Neck ND13</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	
		LP11 09 0627 LP11 09 0644 LP11 09 0645 LP11 14 1190 LP11 19 0933 LP11 19 1022	LP13 09 0222 LP13 09 1335 LP13 09 1336 LP13 09 0280	LP08 14 0641	LP08 14 0513 LP08 14 0595	LP11 14 0544 LP11 14 0545	LP15 14 0548 LP15 14 0562	LP20 09 0342	LP20 09 1405	LP20 09 0297 LP20 09 1223 LP20 09 0440 LP18 09 0864	LP20 09 0795 LP20 09 0796	LP20 09 0873 LP20 09 1222 LP18 09 1306 LP18 09 1310 LP18 09 1307 LP18 09 1311	LP20 09 0801	LP20 09 0802
CTC	PAL								X			X		
CTC	A 200 S	X												
Dani	ALS 39.80	X												
Dani	ALS 86.80	X												
Dani	HS39.50										X			
Dani	HS86.50										X			
Dionex	Autoselect AS50	X												
D-Star	DAS-10													
Dynatech	42 vial tray	X												
Dynatech	60 vial tray													
Dynatech	LC2000													
Finnigan	A 200 S													
Fisons	60 vial tray													
Fisons	AS 800	X												
Fisons	HS500								X					
Fisons	HS800								X					
Fisons	HS850									X				
GBC		X												
Gerstel	MPS2								X			X		
Gilson	(not ASPEC)													
Gilson	201/202					X	X							
Gilson	221/222					X	X							
Gilson	231/401					X								
Gilson	232/402					X								
Gilson	ASPEC					X	X							
Gilson		X												
HTA	HT300A	X												
ICI	LC1600													
Infochroma		X												
KNAUER	Smartline 3800	X				X								
KNAUER	Smartline 3900	X												
Kontron	MSI 660		X											
Kontron		X												
LDC	713-60													
LDC	Marathon	X												



5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Vials

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Rollrand ND8	Rollrand ND8	Rollrand ND8	Rollrand ND8	Rollrand ND8	Gewinde ND8	Gewinde ND8	Kurzgewinde ND9	Gewinde ND10	Rollrand ND11	Rollrand ND11	Rollrand ND11	Rollrand ND11
		<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Screw Neck ND8</i>	<i>Screw Neck ND8</i>	<i>Short Thread ND9</i>	<i>Screw Neck ND10</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>
		LP08 09 0405 LP08 09 0406	LP08 09 0284	LP08 09 0845 LP08 09 0953	LP08 09 0276 LP08 09 0606	LP08 09 0305 LP08 09 0258	LP11 09 0210 LP11 09 0259 LP11 09 0419 LP11 09 0382	LP11 09 0417	LP11 09 0500 LP11 09 0519 LP11 09 0520 LP11 09 0999 LP11 09 0620 LP11 14 1189 LP11 19 0932 LP11 19 1021	LP10 09 0743 LP10 09 1196 LP10 09 1197	LP11 09 0184 LP11 09 0356 LP11 09 0476 LP11 09 0477 LP11 09 0921 LP11 09 0619	LP11 09 0486	LP11 09 0415	LP11 09 0831
LEAP	A200S	X				X								
LEAP	HTS-PAL	X			X	X								
Merck/Hitachi	AS2000								X		X			
Merck/Hitachi	AS4000		X						X		X			
Merck/Hitachi	LaChrom													
Merck/Hitachi					X		X		X					
Perkin Elmer	AI-1		X*						X		X			
Perkin Elmer	AS-100		X											
Perkin Elmer	AS-100B		X											
Perkin Elmer	AS2000													
Perkin Elmer	AS2000B		X											
Perkin Elmer	AS-300		X											
Perkin Elmer	AS8300		X											
Perkin Elmer	Autosystem				X				X		X			
Perkin Elmer	HS6													
Perkin Elmer	HS100													
Perkin Elmer	HS101													
Perkin Elmer	HS40													
Perkin Elmer	Integral4000				X				X		X		X	
Perkin Elmer	ISS-100				X				X		X		X	
Perkin Elmer	ISS-200: 85-, 100-, 145-vial tray				X				X		X		X	
Perkin Elmer	ISS-200: 25-, 85-vial tray				X				X		X			
Perkin Elmer	ISS-200: 145-, 150-vial tray				X				X		X			
Perkin Elmer	LC 600 42 vial tray								X		X			
Perkin Elmer	LC 600 60 vial tray		X											
Perkin Elmer	TurboMatrix HS16													
Perkin Elmer	TurboMatrix HS40													
Perkin Elmer	TurboMatrix HS110													
Pharmacia	LKB 2157-010						X		X		X			
Pharmacia	LKB 2157-020	X				X								X
Quma Elektronik	QHSS-40								X		X			
Sedere									X		X			
SGE	LS-3200	X				X								
Shimadzu	AOC-5000 (Combi PAL)	X			X	X								
Shimadzu	AOC-14/1400				X		X		X		X		X	
Shimadzu	AOC-20						X		X					

* (w. adapter)

5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Flaschen

5.11 Autosampler Compatibility Chart

5. HPLC Vials

Schnapping ND11 <i>Snap Ring ND11</i>	Gewinde ND13 <i>Screw Neck ND13</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)
LP11 09 0627 LP11 09 0644 LP11 09 0645 LP11 14 1190 LP11 19 0933 LP11 19 1022	LP13 09 0222 LP13 09 1335 LP13 09 1336 LP13 09 0280	LP08 14 0641	LP08 14 0513 LP08 14 0595	LP11 14 0544 LP11 14 0545	LP15 14 0548 LP15 14 0562	LP20 09 0342	LP20 09 1405	LP20 09 0297 LP20 09 1223 LP20 09 0440 LP18 09 0864	LP20 09 0795 LP20 09 0796	LP20 09 0873 LP20 09 1222 LP18 09 1306 LP18 09 1310 LP18 09 1307 LP18 09 1311	LP20 09 0801	LP20 09 0802

LEAP	A200S											
LEAP	HTS-PAL						X			X		
Merck/Hitachi	AS2000	X										
Merck/Hitachi	AS4000	X										
Merck/Hitachi	LaChrom		X									
Merck/Hitachi												
Perkin Elmer	AI-1	X										
Perkin Elmer	AS-100											
Perkin Elmer	AS-100B				X							
Perkin Elmer	AS2000				X							
Perkin Elmer	AS2000B				X							
Perkin Elmer	AS-300				X							
Perkin Elmer	AS8300				X							
Perkin Elmer	Autosystem	X	X		X							
Perkin Elmer	HS6					X						
Perkin Elmer	HS100							X				
Perkin Elmer	HS101							X				
Perkin Elmer	HS40							X				
Perkin Elmer	Integral4000	X			X	X						
Perkin Elmer	ISS-100	X			X							
Perkin Elmer	ISS-200: 85-, 100-, 145-vial tray	X			X							
Perkin Elmer	ISS-200: 25-, 85-vial tray			X		X						
Perkin Elmer	ISS-200: 145-, 150-vial tray											
Perkin Elmer	LC 600 42 vial tray	X			X							
Perkin Elmer	LC 600 60 vial tray											
Perkin Elmer	TurboMatrix HS16							X				
Perkin Elmer	TurboMatrix HS40							X				
Perkin Elmer	TurboMatrix HS110							X				
Pharmacia	LKB 2157-010	X										
Pharmacia	LKB 2157-020											
Quma Elektronik	QHSS-40								X			
Sedere		X										
SGE	LS-3200											
Shimadzu	AOC-5000 (Combi PAL)						X			X		
Shimadzu	AOC-14/1400	X	X		X							
Shimadzu	AOC-20											



5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Flaschen

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Rollrand ND8	Rollrand ND8	Rollrand ND8	Rollrand ND8	Rollrand ND8	Gewinde ND8	Gewinde ND8	Kurzgewinde ND9	Gewinde ND10	Rollrand ND11	Rollrand ND11	Rollrand ND11	Rollrand ND11	
		<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Crimp Neck ND8</i>	<i>Screw Neck ND8</i>	<i>Screw Neck ND8</i>	<i>Short Thread ND9</i>	<i>Screw Neck ND10</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	<i>Crimp Neck ND11</i>	
		LP08 09 0405 LP08 09 0406	LP08 09 0284	LP08 09 0845 LP08 09 0953	LP08 09 0276 LP08 09 0606	LP08 09 0305 LP08 09 0258	LP11 09 0210 LP11 09 0259 LP11 09 0419 LP11 09 0382	LP11 09 0417	LP11 09 0500 LP11 09 0519 LP11 09 0520 LP11 09 0999 LP11 09 0620 LP11 14 1189 LP11 19 0932 LP11 19 1021	LP10 09 0743 LP10 09 1196 LP10 09 1197	LP11 09 0184 LP11 09 0356 LP11 09 0476 LP11 09 0477 LP11 09 0921 LP11 09 0619	LP11 09 0486	LP11 09 0415	LP11 09 0831	
Shimadzu	LC-10A				X									X	
Shimadzu	SIL-10A				X		X			X					
Shimadzu	SIL-2AS				X					X					
Shimadzu	SIL-6A				X					X					
Shimadzu	SIL-6B				X		X		X	X	X		X		
Shimadzu	SIL-7A									X					
Shimadzu	SIL-8A									X					
Shimadzu	SIL-9A				X		X		X	X	X		X		
Shimadzu	SIL-HT									X					
Spark	Marathon						X		X		X				
Spark	Midas								X		X				
Spark	Promis						X		X		X				
Spark	SPH 125								X		X				
Spark	Triathlon		X												
Spark					X									X	
Spectra-Physics	8875						X		X		X				
Spectra-Physics	8880						X		X		X				
Spectra-Physics	SpectraSYSTEM								X		X				
Sykam	S 5200								X						
Talbot									X		X				
Tekmar-Dohrmann	7000HT														
Tekmar-Dohrmann	SOLATek72TM								X		X				
ThermoQuest	AS 800						X	X					X	X	
ThermoQuest	HS 500														
ThermoQuest	HS 850														
ThermoQuest	Trace GC 2000 + AS2000						X		X		X				
ThermoQuest	Trace GC 2000 + HS2000														
Tosoh	AS-8010								X		X				
Tosoh	TSK-6080								X		X				
TSP	8875						X		X		X				
Unicam	4247								X		X				
Unicam	4710								X		X				
Unicam	4700 (GC)	X				X									
Unicam	4700 (LC)	X				X	X	X	X						
Unicam	LC-XP						X	X	X		X		X		
Unicam	S4	X				X									

5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Flaschen

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Schnapping ND11 <i>Snap Ring ND11</i>	Gewinde ND13 <i>Screw Neck ND13</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	
		LP11 09 0627 LP11 09 0644 LP11 09 0645 LP11 14 1190 LP11 19 0933 LP11 19 1022	LP13 09 0222 LP13 09 1335 LP13 09 1336 LP13 09 0280	LP08 14 0641	LP08 14 0513 LP08 14 0595	LP11 14 0544 LP11 14 0545	LP15 14 0548 LP15 14 0562	LP20 09 0342	LP20 09 1405	LP20 09 0297 LP20 09 1223 LP20 09 0440 LP18 09 0864	LP20 09 0795 LP20 09 0796	LP20 09 0873 LP20 09 1222 LP18 09 1306 LP18 09 1310 LP18 09 1307 LP18 09 1311	LP20 09 0801	LP20 09 0802
Shimadzu	LC-10A		X											
Shimadzu	SIL-10A					X								
Shimadzu	SIL-2AS		X											
Shimadzu	SIL-6A													
Shimadzu	SIL-6B	X				X								
Shimadzu	SIL-7A					X								
Shimadzu	SIL-8A					X								
Shimadzu	SIL-9A	X				X								
Shimadzu	SIL-HT		X											
Spark	Marathon	X				X								
Spark	Midas	X												
Spark	Promis	X				X								
Spark	SPH 125	X												
Spark	Triathlon					X								
Spark														
Spectra-Physics	8875	X												
Spectra-Physics	8880	X												
Spectra-Physics	SpectraSYSTEM	X												
Sykam	S 5200	X												
Talbot		X												
Tekmar-Dohrmann	7000HT									X				
Tekmar-Dohrmann	SOLATek72TM													
ThermoQuest	AS 800													
ThermoQuest	HS 500							X			X			
ThermoQuest	HS 850													
ThermoQuest	Trace GC 2000 + AS2000	X												
ThermoQuest	Trace GC 2000 + HS2000										X			
Tosoh	AS-8010	X												
Tosoh	TSK-6080	X	X											
TSP	8875	X												
Unicam	4247	X												
Unicam	4710	X												
Unicam	4700 (GC)													
Unicam	4700 (LC)													
Unicam	LC-XP	X	X											
Unicam	S4													



5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Flaschen

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Rollrand ND8 <i>Crimp Neck ND8</i>	Rollrand ND8 <i>Crimp Neck ND8</i>	Rollrand ND8 <i>Crimp Neck ND8</i>	Rollrand ND8 <i>Crimp Neck ND8</i>	Rollrand ND8 <i>Crimp Neck ND8</i>	Gewinde ND8 <i>Screw Neck ND8</i>	Gewinde ND8 <i>Screw Neck ND8</i>	Kurzw Gewinde ND9 <i>Short Thread ND9</i>	Gewinde ND10 <i>Screw Neck ND10</i>	Rollrand ND11 <i>Crimp Neck ND11</i>	Rollrand ND11 <i>Crimp Neck ND11</i>	Rollrand ND11 <i>Crimp Neck ND11</i>	Rollrand ND11 <i>Crimp Neck ND11</i>	
		LP08 09 0405 LP08 09 0406	LP08 09 0284	LP08 09 0845 LP08 09 0953	LP08 09 0276 LP08 09 0606	LP08 09 0305 LP08 09 0258	LP11 09 0210 LP11 09 0259 LP11 09 0419 LP11 09 0382	LP11 09 0417	LP11 09 0500 LP11 09 0519 LP11 09 0520 LP11 09 0999 LP11 09 0620 LP11 14 1189 LP11 19 0932 LP11 19 1021	LP10 09 0743 LP10 09 1196 LP10 09 1197	LP11 09 0184 LP11 09 0356 LP11 09 0476 LP11 09 0477 LP11 09 0921 LP11 09 0619	LP11 09 0486	LP11 09 0415	LP11 09 0831	
Unicam	S8	X				X									
Varian	8035						X	X	X						
Varian	8400														
Varian	8410								X						
Varian	8100						X	X	X		X				
Varian	8200						X	X	X		X				
Varian	9095						X	X	X						
Varian	9100						X	X	X						
Varian	8000 series						X	X	X						
Varian	9090						X	X	X						
Varian	CombiPAL						X	X	X		X				
Varian	Genesis														
Varian	Marathon						X	X	X		X			X	
Varian	Vista						X		X						
Varian (CP)	CP-9020														
Varian (CP)	CP-9025														
Waters	48 position														
Waters	96 position			X											
Waters	2690								X		X				
Waters	Alliance								X	X	X				



5.11 Autosampler Kompatibilitäts-Übersicht

5. HPLC Flaschen

5.11 Autosampler Compatibility Chart

5. HPLC Vials

		Schnapping ND11 <i>Snap Ring ND11</i>	Gewinde ND13 <i>Screw Ring ND13</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Flachboden- gläser <i>Shell Vials</i>	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	Headspace ND20 (ND18)	
		LP11 09 0627 LP11 09 0644 LP11 09 0645 LP11 14 1190 LP11 19 0933 LP11 19 1022	LP13 09 0222 LP13 09 1335 LP13 09 1336 LP13 09 0280	LP08 14 0641	LP08 14 0513 LP08 14 0595	LP11 14 0544 LP11 14 0545	LP15 14 0548 LP15 14 0562	LP20 09 0342	LP20 09 1405	LP20 09 0297 LP20 09 1223 LP20 09 0440 LP18 09 0864	LP20 09 0795 LP20 09 0796	LP20 09 0873 LP20 09 1222 LP18 09 1306 LP18 09 1310 LP18 09 1307 LP18 09 1311	LP20 09 0801	LP20 09 0802
Unicam	S8													
Varian	8035													
Varian	8400													
Varian	8410													
Varian	8100	X									X		X	
Varian	8200	X												
Varian	9095													
Varian	9100													
Varian	8000 series													
Varian	9090													
Varian	CombiPAL	X												
Varian	Genesis							X	X	X	X			
Varian	Marathon	X												
Varian	Vista													
Varian (CP)	CP-9020								X					
Varian (CP)	CP-9025								X					
Waters	48 position		X				X							
Waters	96 position				X									
Waters	2690	X												
Waters	Alliance	X												

Zum Zeitpunkt des Druckes war diese Kompatibilitätsübersicht korrekt. Bitte beachten Sie, dass spätere Änderungen in den Spezifikationen der Hersteller zu abweichenden Kompatibilitäten führen können.

According to our information our compatibility chart was correct when it was printed. Please note that later changes in manufacturers' specifications may lead to different compatibilities.